













Symbol Glossary




Symbol	Language	Title of symbol	Description of symbol	Title & Designation Number of Standard	Symbol Reference Number
General					
	American English	Manufacturer	Indicates the medical device manufacturer	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.1.1.
	Bosnian	Proizvođač			
	British English	Manufacturer			
	Bulgarian	Производител			
	Croatian	Proizvođač			
	Czech	Výrobce			
	Danish	Producent			
	Dutch	Fabrikant			
	Estonian	Tootja			
	Finnish	Valmistaja			
	French	Fabricant			
	German	Hersteller			
	Greek	Κατασκευαστής			
	Hungarian	Gyártó			
	Icelandic	Framleiðandi			
	Italian	Produttore			
	Latvian	Ražotājs			
	Lithuanian	Gamintojas			
	Macedonian	Производител			
	Norwegian	Produsent			
Polish	Producent				
Portuguese	Fabricante				
Romanian	Producător				
Serbian	Proizvođač				
Slovak	Výrobca				
Slovenian	Izdelovalec				
Spanish	Fabricante				
Swedish	Tillverkare				
Turkish	Üretici				
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 2px;">EC</div> <div style="border: 1px solid black; padding: 2px;">REP</div> </div>	American English	Authorized representative in the European Community/European Union	Indicates the authorized representative in the European Community / European Union	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.1.2.
	Bosnian	Ovlašteni predstavnik u Evropskoj zajednici/Evropskoj uniji			
	British English	Authorised representative in the European Community/European Union			
	Bulgarian	Упълномощен представител в Европейската общност / Европейския съюз			
	Croatian	Ovlašteni zastupnik u Evropskoj zajednici / Evropskoj uniji			
	Czech	Zplnomocněný zástupce v Evropském společenství / Evropské unii			
	Danish	Autoriseret repræsentant i det Europæiske Fællesskab/Europæiske Union			
	Dutch	Gemachtigde in de Europese Gemeenschap/Europese Unie			
	Estonian	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus			
	Finnish	Euroopan yhteisön / Euroopan unionin valtuutettu edustaja			
	French	Mandataire établi dans l'Union européenne			
	German	Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union			
	Greek	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Ευρωπαϊκή Ένωση			
	Hungarian	Meghatalmazott képviselő az Európai Közösségekben/Európai Unióban			
	Icelandic	Umboðsaðili innan Evrópusambandsins			
	Italian	Rappresentante autorizzato nella Comunità europea / Unione europea			
	Latvian	Pilvarotais pārstāvis Eiropas Kopienā / Eiropas Savienībā			
	Lithuanian	Igaliojatis atstovas Europos Bendrijoje / Europos Sąjungoje			
	Macedonian	Овластен претставник во Европската заедница/Европската Унија			
	Norwegian	Autorisert representant i EU			
Polish	Autoryzowany przedstawiciel na terenie Wspólnoty Europejskiej / Unii Europejskiej				
Portuguese	Representante autorizado na Comunidade Europeia/União Europeia				
Romanian	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană				
Serbian	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj uniji				
Slovak	Spĺnomocnený zástupca v Európskom spoločenstve alebo Európskej únii				
Slovenian	Pooblaščen zastopnik v Evropski skupnosti/Evropski uniji				
Spanish	Representante autorizado en la Comunidad Europea o la Unión Europea				
Swedish	Auktoriserad representant för Europeiska gemenskapen/Europeiska unionen				
Turkish	Avrupa Topluluğunda / Birliğinde yetkili temsilci				
	American English	Date of manufacture	Indicates the date (YYYY-MM-DD) when the medical device was manufactured	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.1.3.
	Bosnian	Datum proizvodnje			
	British English	Date of manufacture			
	Bulgarian	Дата на производство			
	Croatian	Datum proizvodnje			
	Czech	Datum výroby			
	Danish	Fremstillingsdato			
	Dutch	Productiedatum			
	Estonian	Tootmiskuupäev			
	Finnish	Valmistuspäivämäärä			
	French	Date de fabrication			
	German	Herstellungdatum			
	Greek	Ημερομηνία κατασκευαστή			
	Hungarian	Gyártás dátuma			
	Icelandic	Framleiðsludagsetning			
	Italian	Data di produzione			
	Latvian	Izgatavošanas datums			
	Lithuanian	Pagamavimo data			
	Macedonian	Датум на производство			
	Norwegian	Produksjonsdato			
Polish	Data produkcji				
Portuguese	Data de fabrico				
Romanian	Data fabricării				
Serbian	Datum proizvodnje				
Slovak	Dátum výroby				
Slovenian	Datum izdelave				
Spanish	Fecha de fabricación				
Swedish	Tillverkningsdatum				
Turkish	Üretim tarihi				




	American English	Use-by date	Indicates the date (YYYY-MM-DD) after which the medical device is not to be used	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.1.4.
	Bosnian	Upotrijebite do datuma			
	British English	Use-by date			
	Bulgarian	Срок на годност			
	Croatian	Upotrijebiti do			
	Czech	Datum použitelnosti			
	Danish	Sidste anvendelsesdato			
	Dutch	Te gebruiken tot			
	Estonian	Kasutamise lõpptähtpäev			
	Finnish	Viimeinen käyttöpäivä			
	French	Date de péremption			
	German	Halbbarkeitsdatum			
	Greek	Ημερομηνία λήξης			
	Hungarian	Felhasználható			
	Icelandic	Notist fyrir dags			
	Italian	Data di scadenza			
	Latvian	Izlietot līdz			
	Lithuanian	Tinkamumo naudoti data			
	Macedonian	Краен датум на употреба			
	Norwegian	Utløpsdato			
	Polish	Termin ważności			
	Portuguese	Data de validade			
	Romanian	A se utiliza înainte de			
	Serbian	Datum važenja			
	Slovak	Dátum spotreby			
	Slovenian	Porabiti de			
Spanish	Fecha de caducidad				
Swedish	Använd före datum				
Turkish	Son kullanma tarihi				
	American English	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.1.5.
	Bosnian	Broj serije			
	British English	Batch code			
	Bulgarian	Партиден код			
	Croatian	Šifra serije			
	Czech	Kód šarže			
	Danish	Batchkode			
	Dutch	Lotnummer			
	Estonian	Partii kood			
	Finnish	Eräkoodi			
	French	Code du lot			
	German	Chargennummer			
	Greek	Κωδικός партиδας			
	Hungarian	Gyártási tétel kódja			
	Icelandic	Lotukóði			
	Italian	Lotto			
	Latvian	Partijas kods			
	Lithuanian	Partijos kodas			
	Macedonian	Код на серијата			
	Norwegian	Batchkode			
	Polish	Numer serii			
	Portuguese	Código de lote			
	Romanian	Codul lotului			
	Serbian	Šifra šarže			
	Slovak	Kód dávky			
	Slovenian	Številka serije			
Spanish	Código de lote				
Swedish	Lot-nummer				
Turkish	Parti kodu				
	American English	Catalog number	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.1.6.
	Bosnian	Kataloški broj			
	British English	Catalogue number			
	Bulgarian	Каталожен номер			
	Croatian	Kataloški broj			
	Czech	Katalogové číslo			
	Danish	Katalognummer			
	Dutch	Catalogusnummer			
	Estonian	Katalooginumber			
	Finnish	Luettelonumero			
	French	Référence catalogue			
	German	Katalognummer			
	Greek	Αριθμός καταλόγου			
	Hungarian	Katalógusszám			
	Icelandic	Númer í bæklingi			
	Italian	Numero di catalogo			
	Latvian	Kataloga numurs			
	Lithuanian	Katalogo numeris			
	Macedonian	Каталожки број			
	Norwegian	Katalognummer			
	Polish	Numer katalogowy			
	Portuguese	Número do catálogo			
	Romanian	Număr catalog			
	Serbian	Kataloška oznaka			
	Slovak	Katalogové číslo			
	Slovenian	Kataloška številka			
Spanish	Número de catálogo				
Swedish	Katalognummer				
Turkish	Katalog numarası				




	<table border="1"> <tbody> <tr><td>American English</td><td>Importer</td></tr> <tr><td>Bosnian</td><td>Uvoznik</td></tr> <tr><td>British English</td><td>Importer</td></tr> <tr><td>Bulgarian</td><td>Вносител</td></tr> <tr><td>Croatian</td><td>Uvoznik</td></tr> <tr><td>Czech</td><td>Dovozce</td></tr> <tr><td>Danish</td><td>Importør</td></tr> <tr><td>Dutch</td><td>Importeur</td></tr> <tr><td>Estonian</td><td>Importija</td></tr> <tr><td>Finnish</td><td>Maahantuojaja</td></tr> <tr><td>French</td><td>Importateur</td></tr> <tr><td>German</td><td>Importeur</td></tr> <tr><td>Greek</td><td>Εισαγωγέας</td></tr> <tr><td>Hungarian</td><td>Importőr</td></tr> <tr><td>Icelandic</td><td>Innflytjandi</td></tr> <tr><td>Italian</td><td>Importatore</td></tr> <tr><td>Latvian</td><td>Importētājs</td></tr> <tr><td>Lithuanian</td><td>Importuotojas</td></tr> <tr><td>Macedonian</td><td>Увозник</td></tr> <tr><td>Norwegian</td><td>Importør</td></tr> <tr><td>Polish</td><td>Importer</td></tr> <tr><td>Portuguese</td><td>Importador</td></tr> <tr><td>Romanian</td><td>Importator</td></tr> <tr><td>Serbian</td><td>Uvoznik</td></tr> <tr><td>Slovak</td><td>Dovozca</td></tr> <tr><td>Slovenian</td><td>Uvoznik</td></tr> <tr><td>Spanish</td><td>Importador</td></tr> <tr><td>Swedish</td><td>Importör</td></tr> <tr><td>Turkish</td><td>İthalatçı</td></tr> </tbody> </table>	American English	Importer	Bosnian	Uvoznik	British English	Importer	Bulgarian	Вносител	Croatian	Uvoznik	Czech	Dovozce	Danish	Importør	Dutch	Importeur	Estonian	Importija	Finnish	Maahantuojaja	French	Importateur	German	Importeur	Greek	Εισαγωγέας	Hungarian	Importőr	Icelandic	Innflytjandi	Italian	Importatore	Latvian	Importētājs	Lithuanian	Importuotojas	Macedonian	Увозник	Norwegian	Importør	Polish	Importer	Portuguese	Importador	Romanian	Importator	Serbian	Uvoznik	Slovak	Dovozca	Slovenian	Uvoznik	Spanish	Importador	Swedish	Importör	Turkish	İthalatçı	<p>Indicates the entity importing the medical device into the locale</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.1.8.</p>
American English	Importer																																																													
Bosnian	Uvoznik																																																													
British English	Importer																																																													
Bulgarian	Вносител																																																													
Croatian	Uvoznik																																																													
Czech	Dovozce																																																													
Danish	Importør																																																													
Dutch	Importeur																																																													
Estonian	Importija																																																													
Finnish	Maahantuojaja																																																													
French	Importateur																																																													
German	Importeur																																																													
Greek	Εισαγωγέας																																																													
Hungarian	Importőr																																																													
Icelandic	Innflytjandi																																																													
Italian	Importatore																																																													
Latvian	Importētājs																																																													
Lithuanian	Importuotojas																																																													
Macedonian	Увозник																																																													
Norwegian	Importør																																																													
Polish	Importer																																																													
Portuguese	Importador																																																													
Romanian	Importator																																																													
Serbian	Uvoznik																																																													
Slovak	Dovozca																																																													
Slovenian	Uvoznik																																																													
Spanish	Importador																																																													
Swedish	Importör																																																													
Turkish	İthalatçı																																																													
	<table border="1"> <tbody> <tr><td>American English</td><td>Medical device</td></tr> <tr><td>Bosnian</td><td>Medicinsko sredstvo</td></tr> <tr><td>British English</td><td>Medical device</td></tr> <tr><td>Bulgarian</td><td>Медицинско изделие</td></tr> <tr><td>Croatian</td><td>Medicinski proizvod</td></tr> <tr><td>Czech</td><td>Zdravotnický prostředek</td></tr> <tr><td>Danish</td><td>Medicinsk udstyr</td></tr> <tr><td>Dutch</td><td>Medisch hulpmiddel</td></tr> <tr><td>Estonian</td><td>Meditsiiniseade</td></tr> <tr><td>Finnish</td><td>Lääkinnällinen laite</td></tr> <tr><td>French</td><td>Dispositif médical</td></tr> <tr><td>German</td><td>Medizinprodukt</td></tr> <tr><td>Greek</td><td>Ιατροτεχνολογικό προϊόν</td></tr> <tr><td>Hungarian</td><td>Orvostechnikai eszköz</td></tr> <tr><td>Icelandic</td><td>Lækningateki</td></tr> <tr><td>Italian</td><td>Dispositivo medico</td></tr> <tr><td>Latvian</td><td>Medicīniskā ierīce</td></tr> <tr><td>Lithuanian</td><td>Medicinos priemonė</td></tr> <tr><td>Macedonian</td><td>Медицински уред</td></tr> <tr><td>Norwegian</td><td>Medisinsk utstyr</td></tr> <tr><td>Polish</td><td>Urządzenie medyczne</td></tr> <tr><td>Portuguese</td><td>Dispositivo médico</td></tr> <tr><td>Romanian</td><td>Dispozitiv medical</td></tr> <tr><td>Serbian</td><td>Medicinski proizvod</td></tr> <tr><td>Slovak</td><td>Zdravotnícke pomôcky</td></tr> <tr><td>Slovenian</td><td>Medicinski pripomoček</td></tr> <tr><td>Spanish</td><td>Producto sanitario</td></tr> <tr><td>Swedish</td><td>Medicinteknisk produkt</td></tr> <tr><td>Turkish</td><td>Tıbbi cihaz</td></tr> </tbody> </table>	American English	Medical device	Bosnian	Medicinsko sredstvo	British English	Medical device	Bulgarian	Медицинско изделие	Croatian	Medicinski proizvod	Czech	Zdravotnický prostředek	Danish	Medicinsk udstyr	Dutch	Medisch hulpmiddel	Estonian	Meditsiiniseade	Finnish	Lääkinnällinen laite	French	Dispositif médical	German	Medizinprodukt	Greek	Ιατροτεχνολογικό προϊόν	Hungarian	Orvostechnikai eszköz	Icelandic	Lækningateki	Italian	Dispositivo medico	Latvian	Medicīniskā ierīce	Lithuanian	Medicinos priemonė	Macedonian	Медицински уред	Norwegian	Medisinsk utstyr	Polish	Urządzenie medyczne	Portuguese	Dispositivo médico	Romanian	Dispozitiv medical	Serbian	Medicinski proizvod	Slovak	Zdravotnícke pomôcky	Slovenian	Medicinski pripomoček	Spanish	Producto sanitario	Swedish	Medicinteknisk produkt	Turkish	Tıbbi cihaz	<p>Indicates the item is a medical device</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.7.7.</p>
American English	Medical device																																																													
Bosnian	Medicinsko sredstvo																																																													
British English	Medical device																																																													
Bulgarian	Медицинско изделие																																																													
Croatian	Medicinski proizvod																																																													
Czech	Zdravotnický prostředek																																																													
Danish	Medicinsk udstyr																																																													
Dutch	Medisch hulpmiddel																																																													
Estonian	Meditsiiniseade																																																													
Finnish	Lääkinnällinen laite																																																													
French	Dispositif médical																																																													
German	Medizinprodukt																																																													
Greek	Ιατροτεχνολογικό προϊόν																																																													
Hungarian	Orvostechnikai eszköz																																																													
Icelandic	Lækningateki																																																													
Italian	Dispositivo medico																																																													
Latvian	Medicīniskā ierīce																																																													
Lithuanian	Medicinos priemonė																																																													
Macedonian	Медицински уред																																																													
Norwegian	Medisinsk utstyr																																																													
Polish	Urządzenie medyczne																																																													
Portuguese	Dispositivo médico																																																													
Romanian	Dispozitiv medical																																																													
Serbian	Medicinski proizvod																																																													
Slovak	Zdravotnícke pomôcky																																																													
Slovenian	Medicinski pripomoček																																																													
Spanish	Producto sanitario																																																													
Swedish	Medicinteknisk produkt																																																													
Turkish	Tıbbi cihaz																																																													
	<table border="1"> <tbody> <tr><td>American English</td><td>Indicates a barcode as containing Unique Device Identification</td></tr> <tr><td>Bosnian</td><td>Ukazuje na barkod koji sadrži jedinstvenu identifikacijsku oznaku uređaja</td></tr> <tr><td>British English</td><td>Indicates a barcode as containing Unique Device Identification</td></tr> <tr><td>Bulgarian</td><td>Указва баркод, съдържащ Уникален идентификатор на изделие</td></tr> <tr><td>Croatian</td><td>Označava barkod koji sadrži Jedinstvenu identifikaciju proizvoda</td></tr> <tr><td>Czech</td><td>Označuje čárový kód, který obsahuje identifikaci speciálního zařízení</td></tr> <tr><td>Danish</td><td>Indikerer, at strekkoden indeholder unik udstyrsidentifikation</td></tr> <tr><td>Dutch</td><td>Geeft een streepjescode aan die een unieke hulpmiddelidentificatie bevat</td></tr> <tr><td>Estonian</td><td>Näitab, et vöötкод sisaldab seadme kordumatut tunnust</td></tr> <tr><td>Finnish</td><td>Osoittaa, että viivakoodi sisältää laitteen yksilöidyn tunnusteen</td></tr> <tr><td>French</td><td>Présente un code-barres qui contient le numéro d'identification unique du dispositif</td></tr> <tr><td>German</td><td>Zeigt an, dass ein Barcode eine Unique Device Identification (UDI) enthält</td></tr> <tr><td>Greek</td><td>Δείχνει έναν γραμμωτό κώδικα ο οποίος περιέχει το Μοναδικό αναγνωριστικό συσκευής</td></tr> <tr><td>Hungarian</td><td>Az egyedi eszközzazonosítót tartalmazó vonalkódot jelez</td></tr> <tr><td>Icelandic</td><td>Gefur til kynna strikamerki sem inniheldur auðkenni tækis</td></tr> <tr><td>Italian</td><td>Indikerer en streckkod som innehåller unik produktidentifisering</td></tr> <tr><td>Latvian</td><td>Norāda svitrkodu kā unikālu ierīces identifikāciju</td></tr> <tr><td>Lithuanian</td><td>Nurodo brūkšninį kodą, kuriame yra unikalūs priemonės identifikatoriai</td></tr> <tr><td>Macedonian</td><td>Означува баркод што содржи единствена идентификација на уред</td></tr> <tr><td>Norwegian</td><td>Indikerer en strekkode som inneholder unik enhetsidentifikasjon</td></tr> <tr><td>Polish</td><td>Wskazuje, że kod kreskowy zawiera niepowtarzalny numer identyfikacyjny urządzenia</td></tr> <tr><td>Portuguese</td><td>Indica um código de barras contendo a Identificação Única de Dispositivo</td></tr> <tr><td>Romanian</td><td>Indică un cod de bare drept conținând datele de identificare unică a dispozitivului</td></tr> <tr><td>Serbian</td><td>Označava bar-kod koji sadrži Jedinstvenu identifikaciju proizvoda</td></tr> <tr><td>Slovak</td><td>Označuje čiarový kód, ktorý obsahuje jedinečný identifikátor pomôcky</td></tr> <tr><td>Slovenian</td><td>Nakazuje črtno kodo, ki vsebuje enolični identifikator naprave</td></tr> <tr><td>Spanish</td><td>Indica un código de barras que contiene la Identificación Única del Dispositivo</td></tr> <tr><td>Swedish</td><td>Indikerer en streckkod som innehåller unik produktidentifisering</td></tr> <tr><td>Turkish</td><td>Benzersiz Cihaz Kimliği'ni içeren bir barkodu gösterir</td></tr> </tbody> </table>	American English	Indicates a barcode as containing Unique Device Identification	Bosnian	Ukazuje na barkod koji sadrži jedinstvenu identifikacijsku oznaku uređaja	British English	Indicates a barcode as containing Unique Device Identification	Bulgarian	Указва баркод, съдържащ Уникален идентификатор на изделие	Croatian	Označava barkod koji sadrži Jedinstvenu identifikaciju proizvoda	Czech	Označuje čárový kód, který obsahuje identifikaci speciálního zařízení	Danish	Indikerer, at strekkoden indeholder unik udstyrsidentifikation	Dutch	Geeft een streepjescode aan die een unieke hulpmiddelidentificatie bevat	Estonian	Näitab, et vöötкод sisaldab seadme kordumatut tunnust	Finnish	Osoittaa, että viivakoodi sisältää laitteen yksilöidyn tunnusteen	French	Présente un code-barres qui contient le numéro d'identification unique du dispositif	German	Zeigt an, dass ein Barcode eine Unique Device Identification (UDI) enthält	Greek	Δείχνει έναν γραμμωτό κώδικα ο οποίος περιέχει το Μοναδικό αναγνωριστικό συσκευής	Hungarian	Az egyedi eszközzazonosítót tartalmazó vonalkódot jelez	Icelandic	Gefur til kynna strikamerki sem inniheldur auðkenni tækis	Italian	Indikerer en streckkod som innehåller unik produktidentifisering	Latvian	Norāda svitrkodu kā unikālu ierīces identifikāciju	Lithuanian	Nurodo brūkšninį kodą, kuriame yra unikalūs priemonės identifikatoriai	Macedonian	Означува баркод што содржи единствена идентификација на уред	Norwegian	Indikerer en strekkode som inneholder unik enhetsidentifikasjon	Polish	Wskazuje, że kod kreskowy zawiera niepowtarzalny numer identyfikacyjny urządzenia	Portuguese	Indica um código de barras contendo a Identificação Única de Dispositivo	Romanian	Indică un cod de bare drept conținând datele de identificare unică a dispozitivului	Serbian	Označava bar-kod koji sadrži Jedinstvenu identifikaciju proizvoda	Slovak	Označuje čiarový kód, ktorý obsahuje jedinečný identifikátor pomôcky	Slovenian	Nakazuje črtno kodo, ki vsebuje enolični identifikator naprave	Spanish	Indica un código de barras que contiene la Identificación Única del Dispositivo	Swedish	Indikerer en streckkod som innehåller unik produktidentifisering	Turkish	Benzersiz Cihaz Kimliği'ni içeren bir barkodu gösterir	<p>Indicates a carrier that contains Unique Device Identifier information</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.7.10.</p>
American English	Indicates a barcode as containing Unique Device Identification																																																													
Bosnian	Ukazuje na barkod koji sadrži jedinstvenu identifikacijsku oznaku uređaja																																																													
British English	Indicates a barcode as containing Unique Device Identification																																																													
Bulgarian	Указва баркод, съдържащ Уникален идентификатор на изделие																																																													
Croatian	Označava barkod koji sadrži Jedinstvenu identifikaciju proizvoda																																																													
Czech	Označuje čárový kód, který obsahuje identifikaci speciálního zařízení																																																													
Danish	Indikerer, at strekkoden indeholder unik udstyrsidentifikation																																																													
Dutch	Geeft een streepjescode aan die een unieke hulpmiddelidentificatie bevat																																																													
Estonian	Näitab, et vöötкод sisaldab seadme kordumatut tunnust																																																													
Finnish	Osoittaa, että viivakoodi sisältää laitteen yksilöidyn tunnusteen																																																													
French	Présente un code-barres qui contient le numéro d'identification unique du dispositif																																																													
German	Zeigt an, dass ein Barcode eine Unique Device Identification (UDI) enthält																																																													
Greek	Δείχνει έναν γραμμωτό κώδικα ο οποίος περιέχει το Μοναδικό αναγνωριστικό συσκευής																																																													
Hungarian	Az egyedi eszközzazonosítót tartalmazó vonalkódot jelez																																																													
Icelandic	Gefur til kynna strikamerki sem inniheldur auðkenni tækis																																																													
Italian	Indikerer en streckkod som innehåller unik produktidentifisering																																																													
Latvian	Norāda svitrkodu kā unikālu ierīces identifikāciju																																																													
Lithuanian	Nurodo brūkšninį kodą, kuriame yra unikalūs priemonės identifikatoriai																																																													
Macedonian	Означува баркод што содржи единствена идентификација на уред																																																													
Norwegian	Indikerer en strekkode som inneholder unik enhetsidentifikasjon																																																													
Polish	Wskazuje, że kod kreskowy zawiera niepowtarzalny numer identyfikacyjny urządzenia																																																													
Portuguese	Indica um código de barras contendo a Identificação Única de Dispositivo																																																													
Romanian	Indică un cod de bare drept conținând datele de identificare unică a dispozitivului																																																													
Serbian	Označava bar-kod koji sadrži Jedinstvenu identifikaciju proizvoda																																																													
Slovak	Označuje čiarový kód, ktorý obsahuje jedinečný identifikátor pomôcky																																																													
Slovenian	Nakazuje črtno kodo, ki vsebuje enolični identifikator naprave																																																													
Spanish	Indica un código de barras que contiene la Identificación Única del Dispositivo																																																													
Swedish	Indikerer en streckkod som innehåller unik produktidentifisering																																																													
Turkish	Benzersiz Cihaz Kimliği'ni içeren bir barkodu gösterir																																																													




 	American English	Indicates that the product is in compliance with European legislation for medical devices	<p>A marking by which the manufacturer indicates that a device is in conformity with the applicable requirements set out in the European legislation for medical devices and other applicable Union harmonisation legislation providing for its affixing</p>	<p>European Medical Devices Directive 93/42/EEC of 14 June 1993 (as amended by Directive 2007/47/EC)</p> <p>Medical Device Regulation 2017/745 (repealing of MDD 93/42/EEC and 90/385/EEC)</p>	<p>N/A</p>
	Bosnian	Ukazuje da je proizvod uskladen s evropskim zakonodavstvom o medicinskim sredstvima			
	British English	Indicates that the product is in compliance with European legislation for medical devices			
	Bulgarian	Означава, че продуктът е в съответствие с европейското законодателство за медицински изделия			
	Croatian	Znači da je proizvod sukladan s europskim zakonodavstvom za medicinske proizvode			
	Czech	Označuje, že produkt je v souladu s evropskou legislativou pro zdravotnické prostředky			
	Danish	Indikerer, at produktet overholder den europæiske lovgivning for medicinsk udstyr			
	Dutch	Geeft aan dat het product in overeenstemming is met de Europese wetgeving voor medische hulpmiddelen			
	Estonian	Näitab, et toode vastab meditsiiniseadmeid käsitlevatele Euroopa õigusaktidele			
	Finnish	Osoittaa, että tuote noudattaa Euroopan lääketieteellisiä ja kirurgisia laitteita koskevaa lainsäädäntöä			
	French	Indique que le produit est conforme à la législation européenne concernant les dispositifs médicaux			
	German	Einhaltung der europäischen Bestimmungen für Medizinprodukte			
	Greek	Υποδεικνύει ότι το προϊόν συμμορφώνεται με την ευρωπαϊκή νομοθεσία για τα ιατροτεχνολογικά προϊόντα			
	Hungarian	Jelzi, hogy a termék megfelel az orvostechnikai eszközökre vonatkozó európai jogszabályoknak			
	Icelandic	Gefur til kynna að varan uppfyllir evrópskar kröfur um lækningatæki			
	Italian	Indica che il prodotto è conforme alle norme europee relative ai dispositivi medici			
	Latvian	Norāda izstrādājuma atbilstību Eiropas likumdošanai par medicīniskajām ierīcēm			
	Lithuanian	Nurodo, kad gaminys atitinka Europos įstatymus dėl medicinios priemonių			
	Macedonian	Означува дека производот е во согласност со европското законодавство за медицински уреди			
	Norwegian	Indikerer at produktet er i samsvar med europeisk lovgivning for medisinsk utstyr			
Polish	Opisuje, że produkt jest zgodny z przepisami europejskimi dotyczącymi wyrobów medycznych				
Portuguese	Indica que o produto está em conformidade com a legislação europeia para dispositivos médicos				
Romanian	Indică faptul că produsul este în conformitate cu legislația europeană privind dispozitivele medicale				
Serbian	Označava da je proizvod u skladu sa Evropskim zakonima za medicinske proizvode				
Slovak	Označuje, že produkt je v súlade s európskymi právnymi predpismi o zdravotníckych pomôckach				
Slovenian	Pomeni, da je izdelek skladen z evropsko zakonodajo o medicinskih pripomočkih				
Spanish	Indica que el producto cumple la legislación europea para productos sanitarios				
Swedish	Indikerar att produkten efterlever den europeiska lagstiftningen för medicintekniska produkter				
Turkish	Ürünün tbbi cihazlara ilişkin Avrupa mevzuatına uygun olduğunu işaret eder				
	American English	Global Trade Item Number	<p>Indicates a number used to identify trade items at various packaging levels</p>	<p>N/A</p>	<p>N/A</p>
	Bosnian	Globalni broj trgovačke jedinice			
	British English	Global Trade Item Number			
	Bulgarian	Глобален търговски номер на артикул			
	Croatian	Globalni broj trgovačke jedinice			
	Czech	Globální číslo obchodní položky			
	Danish	GTIN-varenummer			
	Dutch	Wereldwijd handelsartikelnummer			
	Estonian	Ülemaailme kaubaartikli kood			
	Finnish	Kauppanimikkeen maailmanlaajuinen yksilöity numero			
	French	Code article international			
	German	Globale Artikelnummer			
	Greek	Διεθνής Κωδικός Μονάδων Εμπορίας			
	Hungarian	Globális kereskedelmi cikkszám			
	Icelandic	Vörunúmer alþjóðlegrar verslunar			
	Italian	Identificatore di prodotto globale per uso commerciale			
	Latvian	Globālās tirdzniecības izstrādājuma numurs			
	Lithuanian	Visuotinis prekinio vieneto numeris			
	Macedonian	Број на ставка за глобално тргување			
	Norwegian	Globalt handelsvarenummer			
Polish	Globalny numer jednostki handlowej				
Portuguese	Número de Item Comercial Global				
Romanian	Număr global de articol comercial				
Serbian	Broj globalne trgovinske jedinice				
Slovak	Globálne identifikačné číslo obchodnej jednotky				
Slovenian	Globalna trgovinska številka izdelka				
Spanish	Número de artículo comercial internacional				
Swedish	Globalt artikelnummer				
Turkish	Global Ticari Ürün Numarası				
	American English	Compliance with Aerosol Dispenser Directive 75/324/EEC	<p>Indicates that a aerosol dispenser is in compliance with the Aerosol Dispenser Directive</p>	<p>Aerosol Dispenser Directive 75/324/EEC</p>	<p>N/A</p>
	Bosnian	U skladu s Direktivom o aerosolnim raspršivačima 75/324/EEZ			
	British English	Compliance with Aerosol Dispenser Directive 75/324/EEC			
	Bulgarian	Съответствие с Директива 75/324/ЕИО за аерозолните опаковки			
	Croatian	Sukladnost s Direktivom o aerosolnim raspršivačima 75/324/EEZ			
	Czech	Ve shodě se směrnicí pro dávkovače aerosolů 75/324/EHS			
	Danish	Overensstemmelse med aerosoldirektivet 75/324/EØF			
	Dutch	Voldoet aan de aerosolrichtlijn 75/324/EEG			
	Estonian	Vastab aerosoolballoone käsitlevale direktiivile 75/324/EMÜ			
	Finnish	Aerosolisäiliöitä koskevan direktiivin 75/324/ETY mukainen			
	French	Conforme à la directive relative aux générateurs d'aérosols 75/324/CEE			
	German	Entspricht der Richtlinie für Aerosolpackungen 75/324/EWG			
	Greek	Συμμόρφωση με την Οδηγία σχετικά με τις Συσκευές Αερολυμάτων 75/324/ΕΟΚ			
	Hungarian	Megfelelőség az aeroszoladagolókra vonatkozó 75/324/EKG irányelvnek			
	Icelandic	Tilskipun 75/324/EEC um samræmingu laga aðildarráttanna varðandi úðabúna			
	Italian	Conforme alla direttiva europea 75/324/CEE relativa agli aerosol			
	Latvian	Atbilstība Aerosolu izsmidzinātāju direktīvai 75/324/EEK			
	Lithuanian	Atitiktis aerozolių dozatoriaus direktyvai 75/324/EEB			
	Macedonian	Усогласено со Директивата за распрскување аеросоли 75/324/EEC			
	Norwegian	Samsvar med direktivet om aerosolbeholdere 75/324/EEC			
Polish	Zgodność z Dyrektywą dotyczącą wyrobów aerozolowych 75/324/EWG				
Portuguese	Conformidade com a Directiva de Difusores de Aerosol 75/324/EEC				
Romanian	Conformitate cu Directiva 75/324/CEE privind generatoarele de aerosoli				
Serbian	Uskladenost sa Direktivom o aerosolnim raspršivačima 75/324/EEC				
Slovak	Súlad so smernicou 75/324/EHS týkajúcou sa aerosolových rozprašovačov				
Slovenian	Skladno z Direktivo o aerosolnih razpršilnikih 75/324/EGS				
Spanish	Cumple la Directiva sobre generadores de aerosoles 75/324/CEE				
Swedish	Uppfyller kraven i aerosoldirektivet 75/324/EEG				
Turkish	75/324 / EEC Aerosol Dağıtıcı Direktifi ile Uygunluk				




Safe use


	<table border="1"> <tbody> <tr><td>American English</td><td>Do not re-use</td></tr> <tr><td>Bosnian</td><td>Nije za višekratnu upotrebu</td></tr> <tr><td>British English</td><td>Do not re-use</td></tr> <tr><td>Bulgarian</td><td>Да не се използва повторно</td></tr> <tr><td>Croatian</td><td>Nemojte ponovo upotrebljavati</td></tr> <tr><td>Czech</td><td>Určeno na jedno použití</td></tr> <tr><td>Danish</td><td>Må ikke genbruges</td></tr> <tr><td>Dutch</td><td>Niet voor hergebruik</td></tr> <tr><td>Estonian</td><td>Korduvalt mitte kasutada</td></tr> <tr><td>Finnish</td><td>Ei saa käyttää uudelleen</td></tr> <tr><td>French</td><td>Ne pas réutiliser</td></tr> <tr><td>German</td><td>Nur für den einmaligen Gebrauch</td></tr> <tr><td>Greek</td><td>Μην το επαναχρησιμοποιείτε</td></tr> <tr><td>Hungarian</td><td>Újrafelhasználása tilos</td></tr> <tr><td>Icelandic</td><td>Notaðu aðeins einu sinni</td></tr> <tr><td>Italian</td><td>Non riutilizzare</td></tr> <tr><td>Latvian</td><td>Nelietot atkārtoti</td></tr> <tr><td>Lithuanian</td><td>Nenaudoti pakartotinai</td></tr> <tr><td>Macedonian</td><td>Да не се употребува повторно</td></tr> <tr><td>Norwegian</td><td>Ikke til gjenbruk</td></tr> <tr><td>Polish</td><td>Do jednorazowego użytku</td></tr> <tr><td>Portuguese</td><td>Não reutilizar</td></tr> <tr><td>Romanian</td><td>A nu se reutiliza</td></tr> <tr><td>Serbian</td><td>Ne koristiti ponovo</td></tr> <tr><td>Slovak</td><td>Nepoužívajte opakovane</td></tr> <tr><td>Slovenian</td><td>Ni za ponovno uporabo</td></tr> <tr><td>Spanish</td><td>No reutilizar</td></tr> <tr><td>Swedish</td><td>Endast för engångsbruk</td></tr> <tr><td>Turkish</td><td>Tekrar kullanmayın</td></tr> </tbody> </table>	American English	Do not re-use	Bosnian	Nije za višekratnu upotrebu	British English	Do not re-use	Bulgarian	Да не се използва повторно	Croatian	Nemojte ponovo upotrebljavati	Czech	Určeno na jedno použití	Danish	Må ikke genbruges	Dutch	Niet voor hergebruik	Estonian	Korduvalt mitte kasutada	Finnish	Ei saa käyttää uudelleen	French	Ne pas réutiliser	German	Nur für den einmaligen Gebrauch	Greek	Μην το επαναχρησιμοποιείτε	Hungarian	Újrafelhasználása tilos	Icelandic	Notaðu aðeins einu sinni	Italian	Non riutilizzare	Latvian	Nelietot atkārtoti	Lithuanian	Nenaudoti pakartotinai	Macedonian	Да не се употребува повторно	Norwegian	Ikke til gjenbruk	Polish	Do jednorazowego użytku	Portuguese	Não reutilizar	Romanian	A nu se reutiliza	Serbian	Ne koristiti ponovo	Slovak	Nepoužívajte opakovane	Slovenian	Ni za ponovno uporabo	Spanish	No reutilizar	Swedish	Endast för engångsbruk	Turkish	Tekrar kullanmayın	<p>Indicates a medical device that is intended for one single use only</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.4.2.</p>
American English	Do not re-use																																																													
Bosnian	Nije za višekratnu upotrebu																																																													
British English	Do not re-use																																																													
Bulgarian	Да не се използва повторно																																																													
Croatian	Nemojte ponovo upotrebljavati																																																													
Czech	Určeno na jedno použití																																																													
Danish	Må ikke genbruges																																																													
Dutch	Niet voor hergebruik																																																													
Estonian	Korduvalt mitte kasutada																																																													
Finnish	Ei saa käyttää uudelleen																																																													
French	Ne pas réutiliser																																																													
German	Nur für den einmaligen Gebrauch																																																													
Greek	Μην το επαναχρησιμοποιείτε																																																													
Hungarian	Újrafelhasználása tilos																																																													
Icelandic	Notaðu aðeins einu sinni																																																													
Italian	Non riutilizzare																																																													
Latvian	Nelietot atkārtoti																																																													
Lithuanian	Nenaudoti pakartotinai																																																													
Macedonian	Да не се употребува повторно																																																													
Norwegian	Ikke til gjenbruk																																																													
Polish	Do jednorazowego użytku																																																													
Portuguese	Não reutilizar																																																													
Romanian	A nu se reutiliza																																																													
Serbian	Ne koristiti ponovo																																																													
Slovak	Nepoužívajte opakovane																																																													
Slovenian	Ni za ponovno uporabo																																																													
Spanish	No reutilizar																																																													
Swedish	Endast för engångsbruk																																																													
Turkish	Tekrar kullanmayın																																																													
	<table border="1"> <tbody> <tr><td>American English</td><td>Consult instructions for use</td></tr> <tr><td>Bosnian</td><td>Pogledajte upite za upotrebu</td></tr> <tr><td>British English</td><td>Consult instructions for use</td></tr> <tr><td>Bulgarian</td><td>Направете справка в инструкциите за употреба</td></tr> <tr><td>Croatian</td><td>Pročitajte upute za upotrebu</td></tr> <tr><td>Czech</td><td>Prečtěte si návod k použití</td></tr> <tr><td>Danish</td><td>Se brugervejledning</td></tr> <tr><td>Dutch</td><td>Gebruiksaanwijzing raadplegen</td></tr> <tr><td>Estonian</td><td>Tutvuge kasutusjuhendiga</td></tr> <tr><td>Finnish</td><td>Lue käyttöohje</td></tr> <tr><td>French</td><td>Consultez le mode d'emploi</td></tr> <tr><td>German</td><td>Gebrauchsanweisung befolgen</td></tr> <tr><td>Greek</td><td>Συμβουλευτείτε τις οδηγίες χρήσης</td></tr> <tr><td>Hungarian</td><td>Olvassa el a használati utasítást</td></tr> <tr><td>Icelandic</td><td>Lestu leiðbeiningar fyrir notkun</td></tr> <tr><td>Italian</td><td>Vedere le istruzioni per l'uso</td></tr> <tr><td>Latvian</td><td>Skatīt lietošanas pamācību</td></tr> <tr><td>Lithuanian</td><td>Žr. naudojimo instrukciją</td></tr> <tr><td>Macedonian</td><td>Видете ги инструкциите за употреба</td></tr> <tr><td>Norwegian</td><td>Se bruksanvisningen</td></tr> <tr><td>Polish</td><td>Przeczytaj instrukcję użycia</td></tr> <tr><td>Portuguese</td><td>Consulte as instruções de utilização</td></tr> <tr><td>Romanian</td><td>Consultați instrucțiunile de utilizare</td></tr> <tr><td>Serbian</td><td>Pogledati uputstvo za upotrebu</td></tr> <tr><td>Slovak</td><td>Konzultujte návod na použitie</td></tr> <tr><td>Slovenian</td><td>Glejte navodila za uporabo</td></tr> <tr><td>Spanish</td><td>Consultar las instrucciones de uso</td></tr> <tr><td>Swedish</td><td>Ta del av bruksanvisningen</td></tr> <tr><td>Turkish</td><td>Kullanma talimatına başvurun</td></tr> </tbody> </table>	American English	Consult instructions for use	Bosnian	Pogledajte upite za upotrebu	British English	Consult instructions for use	Bulgarian	Направете справка в инструкциите за употреба	Croatian	Pročitajte upute za upotrebu	Czech	Prečtěte si návod k použití	Danish	Se brugervejledning	Dutch	Gebruiksaanwijzing raadplegen	Estonian	Tutvuge kasutusjuhendiga	Finnish	Lue käyttöohje	French	Consultez le mode d'emploi	German	Gebrauchsanweisung befolgen	Greek	Συμβουλευτείτε τις οδηγίες χρήσης	Hungarian	Olvassa el a használati utasítást	Icelandic	Lestu leiðbeiningar fyrir notkun	Italian	Vedere le istruzioni per l'uso	Latvian	Skatīt lietošanas pamācību	Lithuanian	Žr. naudojimo instrukciją	Macedonian	Видете ги инструкциите за употреба	Norwegian	Se bruksanvisningen	Polish	Przeczytaj instrukcję użycia	Portuguese	Consulte as instruções de utilização	Romanian	Consultați instrucțiunile de utilizare	Serbian	Pogledati uputstvo za upotrebu	Slovak	Konzultujte návod na použitie	Slovenian	Glejte navodila za uporabo	Spanish	Consultar las instrucciones de uso	Swedish	Ta del av bruksanvisningen	Turkish	Kullanma talimatına başvurun	<p>Indicates the need for the user to consult the instructions for use</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.4.3.</p>
American English	Consult instructions for use																																																													
Bosnian	Pogledajte upite za upotrebu																																																													
British English	Consult instructions for use																																																													
Bulgarian	Направете справка в инструкциите за употреба																																																													
Croatian	Pročitajte upute za upotrebu																																																													
Czech	Prečtěte si návod k použití																																																													
Danish	Se brugervejledning																																																													
Dutch	Gebruiksaanwijzing raadplegen																																																													
Estonian	Tutvuge kasutusjuhendiga																																																													
Finnish	Lue käyttöohje																																																													
French	Consultez le mode d'emploi																																																													
German	Gebrauchsanweisung befolgen																																																													
Greek	Συμβουλευτείτε τις οδηγίες χρήσης																																																													
Hungarian	Olvassa el a használati utasítást																																																													
Icelandic	Lestu leiðbeiningar fyrir notkun																																																													
Italian	Vedere le istruzioni per l'uso																																																													
Latvian	Skatīt lietošanas pamācību																																																													
Lithuanian	Žr. naudojimo instrukciją																																																													
Macedonian	Видете ги инструкциите за употреба																																																													
Norwegian	Se bruksanvisningen																																																													
Polish	Przeczytaj instrukcję użycia																																																													
Portuguese	Consulte as instruções de utilização																																																													
Romanian	Consultați instrucțiunile de utilizare																																																													
Serbian	Pogledati uputstvo za upotrebu																																																													
Slovak	Konzultujte návod na použitie																																																													
Slovenian	Glejte navodila za uporabo																																																													
Spanish	Consultar las instrucciones de uso																																																													
Swedish	Ta del av bruksanvisningen																																																													
Turkish	Kullanma talimatına başvurun																																																													
	<table border="1"> <tbody> <tr><td>American English</td><td>Caution</td></tr> <tr><td>Bosnian</td><td>Oprez</td></tr> <tr><td>British English</td><td>Caution</td></tr> <tr><td>Bulgarian</td><td>Внимание</td></tr> <tr><td>Croatian</td><td>Oprez</td></tr> <tr><td>Czech</td><td>Upozornění</td></tr> <tr><td>Danish</td><td>Forsigtig</td></tr> <tr><td>Dutch</td><td>Let op</td></tr> <tr><td>Estonian</td><td>Ettevaatus</td></tr> <tr><td>Finnish</td><td>Varoitus</td></tr> <tr><td>French</td><td>Avertissement</td></tr> <tr><td>German</td><td>Achtung</td></tr> <tr><td>Greek</td><td>Προφύλαξη</td></tr> <tr><td>Hungarian</td><td>Vigyázat</td></tr> <tr><td>Icelandic</td><td>Varðbarráðstafanir</td></tr> <tr><td>Italian</td><td>Precauzioni</td></tr> <tr><td>Latvian</td><td>Uzmanību</td></tr> <tr><td>Lithuanian</td><td>Dėmesio</td></tr> <tr><td>Macedonian</td><td>Предупредување</td></tr> <tr><td>Norwegian</td><td>Advarsel</td></tr> <tr><td>Polish</td><td>Uwaga</td></tr> <tr><td>Portuguese</td><td>Cuidado</td></tr> <tr><td>Romanian</td><td>Atenție</td></tr> <tr><td>Serbian</td><td>Oprez</td></tr> <tr><td>Slovak</td><td>Upozornenie</td></tr> <tr><td>Slovenian</td><td>Svarilo</td></tr> <tr><td>Spanish</td><td>Precaución</td></tr> <tr><td>Swedish</td><td>Försiktighetsåtgärder</td></tr> <tr><td>Turkish</td><td>Önlem</td></tr> </tbody> </table>	American English	Caution	Bosnian	Oprez	British English	Caution	Bulgarian	Внимание	Croatian	Oprez	Czech	Upozornění	Danish	Forsigtig	Dutch	Let op	Estonian	Ettevaatus	Finnish	Varoitus	French	Avertissement	German	Achtung	Greek	Προφύλαξη	Hungarian	Vigyázat	Icelandic	Varðbarráðstafanir	Italian	Precauzioni	Latvian	Uzmanību	Lithuanian	Dėmesio	Macedonian	Предупредување	Norwegian	Advarsel	Polish	Uwaga	Portuguese	Cuidado	Romanian	Atenție	Serbian	Oprez	Slovak	Upozornenie	Slovenian	Svarilo	Spanish	Precaución	Swedish	Försiktighetsåtgärder	Turkish	Önlem	<p>Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.4.4.</p>
American English	Caution																																																													
Bosnian	Oprez																																																													
British English	Caution																																																													
Bulgarian	Внимание																																																													
Croatian	Oprez																																																													
Czech	Upozornění																																																													
Danish	Forsigtig																																																													
Dutch	Let op																																																													
Estonian	Ettevaatus																																																													
Finnish	Varoitus																																																													
French	Avertissement																																																													
German	Achtung																																																													
Greek	Προφύλαξη																																																													
Hungarian	Vigyázat																																																													
Icelandic	Varðbarráðstafanir																																																													
Italian	Precauzioni																																																													
Latvian	Uzmanību																																																													
Lithuanian	Dėmesio																																																													
Macedonian	Предупредување																																																													
Norwegian	Advarsel																																																													
Polish	Uwaga																																																													
Portuguese	Cuidado																																																													
Romanian	Atenție																																																													
Serbian	Oprez																																																													
Slovak	Upozornenie																																																													
Slovenian	Svarilo																																																													
Spanish	Precaución																																																													
Swedish	Försiktighetsåtgärder																																																													
Turkish	Önlem																																																													

	American English	Contains or presence of natural rubber latex	Indicates the presence of dry natural rubber or natural rubber latex as a material of construction within the medical device or the packaging of a medical device.	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.4.5.
	Bosnian	Sadrži lateks od prirodne gume ili njegove tragove			
	British English	Contains or presence of natural rubber latex			
	Bulgarian	Съдържа или наличие на естествен каучуков латекс			
	Croatian	Sadrži ili je prisutan prirodan gumeni lateks			
	Czech	Obsah nebo přítomnost přírodního kaučukového latexu			
	Danish	Denne produkt indeholder eller kan indeholde naturgummlatex			
	Dutch	Bevat latex, of sporen van latex op basis van natuurlijk rubber			
	Estonian	Loodusliku kummilateksi sisaldus või olemasolu			
	Finnish	Sisältää tai on ollut kosketuksissa luonnonkumilateksin kanssa			
	French	Contient du latex de caoutchouc naturel			
	German	Enthält Naturkautschuklatex			
	Greek	Περιέχει καόλυμα ή παρουσία φυσικό λάτεξ			
	Hungarian	Természetes gumilatexet tartalmaz, vagy az van benne jelen			
	Icelandic	Inniheldur náttúrulegt gúmmilátex, eða það er til staðar			
	Italian	Contenuto o presenza di lattice di gomma naturale			
	Latvian	Satur dabiskā kaučuka lateksu			
	Lithuanian	Sudėtyje yra natūraliojo kaučiuko latekso			
	Macedonian	Содржи или има присуство на природен латекс			
	Norwegian	Inneholder eller kan ha spor av naturgummlateks			
	Polish	Zawiera lateks naturalny			
	Portuguese	Contém ou apresenta látex de borracha natural			
	Romanian	Conținut sau prezentă de latex de cauciuc natural			
	Serbian	Sadrži prirodan lateks			
	Slovak	Obsahuje prírodný latex alebo označuje jeho prítomnosť			
Slovenian	Vsebnost ali prisotnost lateksa iz naravnega kavčuka				
Spanish	Contiene goma de látex natural				
Swedish	Innehåller eller har förekomst av naturgummlatex				
Turkish	Doğal kauçuk lateks içerir veya bulunur				
	American English	Contains a medicinal substance	Indicates a medical device that contains or incorporates a medicinal substance	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.4.7.
	Bosnian	Sadrži ljekovitu supstancu			
	British English	Contains a medicinal substance			
	Bulgarian	Съдържа лекарствено вещество			
	Croatian	Sadrži ljekovitu tvar			
	Czech	Obsahuje léčivou látku			
	Danish	Indeholder et medicinsk stof			
	Dutch	Bevat een medicinale substantie			
	Estonian	Sisaldab meditsiinilist ainet			
	Finnish	Sisältää lääkeainetta			
	French	Contient une substance médicamenteuse			
	German	Enthält eine Arzneimittelsubstanz			
	Greek	Περιέχει ιατρική ουσία			
	Hungarian	Orvostechnikai anyagot tartalmaz			
	Icelandic	Inniheldur lyf			
	Italian	Contiene una sostanza medicinale			
	Latvian	Satur medicīnisku vielu			
	Lithuanian	Sudėtyje yra vaistinės medžiagos			
	Macedonian	Содржи медицинска супстанција			
	Norwegian	Inneholder et medisinsk stoff			
	Polish	Zawiera substancję leczniczą			
	Portuguese	Contém uma substância medicinal			
	Romanian	Contine o substantă medicamentoasă			
	Serbian	Sadrži medicinsku supstancu			
	Slovak	Obsahuje liečivý prípravok			
Slovenian	Vsebuje zdravilno snov				
Spanish	Contiene un principio activo				
Swedish	Innehåller en medicinsk substans				
Turkish	Tibbi madde içerir				
	American English	Contains biological material of animal origin	Indicates a medical device that contains biological tissue, cells, or their derivatives, of animal origin	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.4.8.
	Bosnian	Sadrži biološki materijal životinjskog porijekla			
	British English	Contains biological material of animal origin			
	Bulgarian	Съдържа биологичен материал от животински произход			
	Croatian	Sadrži biološki materijal životinjskog podrijetla			
	Czech	Obsahuje biologický materiál zvířecího původu			
	Danish	Indeholder biologisk materiale af animalsk oprindelse			
	Dutch	Bevat biologisch materiaal van dierlijke oorsprong			
	Estonian	Sisaldab loomset päritolu bioloogilist materjali			
	Finnish	Sisältää elämistä peräisin olevaa biologista ainetta			
	French	Contient un matériau biologique d'origine animale			
	German	Enthält biologisches Material tierischen Ursprungs			
	Greek	Περιέχει βιολογικό υλικό ζωικής προέλευσης			
	Hungarian	Állati eredetű biológiai anyagot tartalmaz			
	Icelandic	Inniheldur lifrænt efni sem kemur úr dýraríkinu			
	Italian	Contiene materiale biologico di origine animale			
	Latvian	Satur dzīvnieku izcelsmes bioloģisko materiālu			
	Lithuanian	Sudėtyje yra gyvūninės kilmės biologinės medžiagos			
	Macedonian	Содржи биолошки материјал од животинско потекло			
	Norwegian	Inneholder biologisk materiale av animalsk opprinnelse			
	Polish	Zawiera materiał biologiczny pochodzenia zwierzęcego			
	Portuguese	Contém material biológico de origem animal			
	Romanian	Conține material biologic de origine animală			
	Serbian	Sadrži biološki materijal životinjskog porekla			
	Slovak	Obsahuje biologický materiál živočíšneho pôvodu			
Slovenian	Vsebuje biološki material živalskega izvora				
Spanish	Contiene material biológico de origen animal				
Swedish	Innehåller biologiskt material av animaliskt ursprung				
Turkish	Hayvan kaynaklı biyolojik materyal içerir				


	American English	Contains hazardous substance	<p>Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR), or substances with endocrine disrupting properties</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.4.10.</p>
	Bosnian	Sadrži opasnu tvar			
	British English	Contains hazardous substance			
	Bulgarian	Съдържа опасно вещество			
	Croatian	Sadrži opasnu tvar			
	Czech	Obsahuje nebezpečnou látku			
	Danish	Indeholder farligt stof			
	Dutch	Bevat een gevaarlijke stof			
	Estonian	Sisaldab ohtlikku ainet			
	Finnish	Sisältää vaarallista ainetta			
	French	Contient une substance dangereuse			
	German	Enthält Gefahrstoff			
	Greek	Περιέχει επικίνδυνα ουσία			
	Hungarian	Veszélyes anyagot tartalmaz			
	Icelandic	Inniheldur hættulegt efni			
	Italian	Contiene una sostanza pericolosa			
	Latvian	Satur bīstamu vielu			
	Lithuanian	Sudėtyje yra pavojingos medžiagos			
	Macedonian	Содржи опасна супстанција			
	Norwegian	Inneholder farlig stoff			
	Polish	Zawiera substancję niebezpieczną			
	Portuguese	Contém substância perigosa			
	Romanian	Contine substanțe periculoase			
Serbian	Sadrži opasne supstance				
Slovak	Obsahuje nebezpečné látky				
Slovenian	Vsebuje nevarno snov				
Spanish	Contiene una sustancia peligrosa				
Swedish	Innehåller farliga ämnen				
Turkish	Tehlikeli madde içerir				
	American English	Single patient multiple use	<p>Indicates a medical device that may be used multiple times (multiple procedures) on a single patient</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.4.12.</p>
	Bosnian	Višestruka upotreba na jednom pacijentu			
	British English	Single patient multiple use			
	Bulgarian	Многokratна употреба при един пациент			
	Croatian	Za višestruku upotrebu na jednom pacijentu			
	Czech	K vícečetnému použití u jednoho pacienta			
	Danish	Flergangsbrug til én enkelt patient			
	Dutch	Meervoudig gebruik voor één patiënt			
	Estonian	Ühel patsiendil korduskasutatav			
	Finnish	Moninkertaiseen käyttöön yhdellä potilaalla			
	French	Peut être utilisé plusieurs fois sur le même patient			
	German	Einzelner Patient – mehrfach verwendbar			
	Greek	Για πολλαπλές χρήσεις σε έναν ασθενή			
	Hungarian	Egyetlen beteg esetében többször újrahasználható			
	Icelandic	Má nota mörgum sinnum hjá einum og sama sjúklingi			
	Italian	Uso multiplo per un solo paziente			
	Latvian	Vairākkārtējai lietošanai vienam pacientam			
	Lithuanian	Vienas pacientas, daugkartinis naudojimas			
	Macedonian	За еден пациент, повеќекратна употреба			
	Norwegian	Til flergangsbruk på kun én pasient			
	Polish	Do wielokrotnego stosowania u jednego pacjenta			
	Portuguese	Único paciente múltiplas utilizações			
	Romanian	Utilizare multiplă la un singur pacient			
Serbian	Višestruka upotreba na jednom pacijentu				
Slovak	Učlenené na viacnásobné použitie pre jedného pacienta				
Slovenian	Večkratna uporaba za enega pacienta				
Spanish	Varios usos en un solo paciente				
Swedish	Avsedd att användas flera gånger på en patient				
Turkish	Tek hastada birden fazla kullanılm				
	American English	Not made with natural rubber latex	<p>Indicates that the medical device is not made with natural rubber latex</p>	<p>N/A</p>	<p>N/A</p>
	Bosnian	Nije proizvedeno od prirodne lateks gume			
	British English	Not made with natural rubber latex			
	Bulgarian	Не е произведено с естествен каучуков латекс			
	Croatian	Nije izrađeno od prirodnog gumenog lateksa			
	Czech	Neobsahuje latex			
	Danish	Ikke fremstillet af naturlig gummilatex			
	Dutch	Niet gemaakt van latex van natuurlijk rubber			
	Estonian	Valmistamisel ei ole kasutatud looduslikku kummilateksit			
	Finnish	Valmistuksessa ei ole käytetty luonnonkumilateksia			
	French	Fabriqu� sans latex de caoutchouc naturel			
	German	Ohne Naturkautschuklatex hergestellt			
	Greek	Δεν είναι κατασκευασμένο από φυσικό λάτεξ			
	Hungarian	Nem tartalmaz term�szetes gumilatexet			
	Icelandic	Ekki framleitt �r n�tt�rulegu g�mmilatexi			
	Italian	Non contiene lattice di gomma naturale			
	Latvian	Nav izgatavots no dabisk� kau�uka lateksa			
	Lithuanian	Pagaminta be nat�raliojo kau�uko latekso			
	Macedonian	He e napraveno co prirodni lareks			
	Norwegian	Ikke produsert av naturlig gummilateks			
	Polish	Wyprodukowano bez użycia naturalnej gumy lateksowej			
	Portuguese	N�o fabricado com l�tex de borracha natural			
	Romanian	Nu este fabricat cu l�tex de cauciuc natural			
Serbian	Nije napravljeno od prirodnog lateksa				
Slovak	Neobsahuje pr�rodn� latex				
Slovenian	Ni izdelano iz lateksa iz naravnega kav�uka				
Spanish	Este producto no est� fabricado con goma de l�tex natural				
Swedish	Inte tillverkad av naturgummilatex				
Turkish	Dođal kau�uk lateks ile �retilememiştir				


	American English	MR safe	<p>Indicates an item that poses no known hazards resulting from exposure to any MR environment. MR Safe items are composed of materials that are electrically nonconductive, nonmetallic, and nonmagnetic</p>	<p>ASTM F 2503 - 20 Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment</p>	<p>3.1.13</p>
	Bosnian	Sigurno za snimanje MR-om			
	British English	MR safe			
	Bulgarian	Безопасно в МР среда			
	Croatian	Sigurno za snimanje MR-om			
	Czech	Bezpečný při vyšetření MR			
	Danish	MR-sikker			
	Dutch	MRI-veilig			
	Estonian	MR-ohutu			
	Finnish	Turvallinen MRI-laitteen kanssa			
	French	Compatible IRM			
	German	MR-sicher			
	Greek	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας			
	Hungarian	MR-biztonságos			
	Icelandic	Öruggt í segulómskoðun (MR)			
	Italian	Sicuro per la RM			
	Latvian	MR drošs			
	Lithuanian	Saugi MR aplinkoje			
	Macedonian	Безбедно за МР			
	Norwegian	MR-sikker			
	Polish	Produkt może przebywać w pobliżu rezonansu magnetycznego			
	Portuguese	Seguro para RM			
	Romanian	Sigur în mediul RM			
Serbian	Bezbedno u okruženju magnetne rezonance				
Slovak	Bezpečné v prostredí MR				
Slovenian	Varno pri MR				
Spanish	Totalmente compatible con RM				
Swedish	MR-säker				
Turkish	MR Güvenli				
	American English	MR conditional	<p>Indicates an item with demonstrated safety in the MR environment within defined conditions for the static magnetic field, the time-varying gradient magnetic fields and the radiofrequency fields</p>	<p>ASTM F 2503 - 20 Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment</p>	<p>3.1.11</p>
	Bosnian	Uslovno sigurno za snimanje MR-om			
	British English	MR conditional			
	Bulgarian	Безопасно в МР среда при определени условия			
	Croatian	Uvjetno sigurno za snimanje MR-om			
	Czech	Podmíněně bezpečný při vyšetření MR			
	Danish	MR-betinget			
	Dutch	MRI-voorwaardelijk			
	Estonian	MR-tingimuslik			
	Finnish	Tietyn ehdoin turvallinen MRI-laitteen kanssa			
	French	Compatible IRM sous conditions			
	German	Bedingt MR-sicher			
	Greek	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας υπό προϋποθέσεις			
	Hungarian	MR-feltételes			
	Icelandic	Þarf að athuga í segulómskoðun (MR)			
	Italian	A compatibilità RM condizionata			
	Latvian	MR ar nosacījumiem			
	Lithuanian	Salyginai saugi MR aplinkoje			
	Macedonian	Условно за МР			
	Norwegian	MR-betinget			
	Polish	Produkt warunkowo może przebywać w pobliżu rezonansu magnetycznego			
	Portuguese	Condicional para RM			
	Romanian	Compatibilitate condiționată în mediul RM			
Serbian	Uslovno bezbedno u okruženju magnetne rezonance				
Slovak	Podmienečne použiteľné v prostredí MR				
Slovenian	Pogojno varno pri MR				
Spanish	Compatibilidad condicional con RM				
Swedish	MR Conditional förutsatt att vissa villkor är uppfyllda				
Turkish	MR Koşullu				
	American English	MR unsafe	<p>Indicates an item which poses unacceptable risks to the patient, medical staff or other persons within the MR environment</p>	<p>ASTM F 2503 - 20 Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment</p>	<p>3.1.14</p>
	Bosnian	Nije sigurno za snimanje MR-om			
	British English	MR unsafe			
	Bulgarian	Не е безопасно в МР среда			
	Croatian	Nije sigurno za snimanje MR-om			
	Czech	Výrobek není bezpečný v prostředí MR			
	Danish	Ikke MR-sikker			
	Dutch	MRI-onveilig			
	Estonian	MR-ohulik			
	Finnish	Epäturvallinen MRI-laitteen kanssa			
	French	Non compatible IRM			
	German	MR-unsicher			
	Greek	Μη ασφαλές σε περιβάλλον μαγνητικής τομογραφίας			
	Hungarian	Nem MR-biztonságos			
	Icelandic	Ekki öruggt í segulómskoðun (MR)			
	Italian	Non compatibile con la RM			
	Latvian	MR nedrošs			
	Lithuanian	Nesaugi MR aplinkoje			
	Macedonian	Не е безбедно за МР			
	Norwegian	MR-usikker			
	Polish	Produkt nie może przebywać w pobliżu rezonansu magnetycznego			
	Portuguese	Não seguro para RM			
	Romanian	Nu este sigur în mediul RM			
Serbian	Nije bezbedno u okruženju magnetne rezonance				
Slovak	Nebezpečné v prostredí MR				
Slovenian	Ni varno pri MR				
Spanish	Totalmente incompatible con RM				
Swedish	Ej MR-säker				
Turkish	MR Güvenli Değil				




	American English	Highly flammable	Indicates a medical device that is highly flammable	United Nations Economic Commission for Europe (UNECE), Globally Harmonized System of Classification and Labelling of Chemicals (GHS)	N/A
	Bosnian	Lako zapaljivo			
	British English	Highly flammable			
	Bulgarian	Силно запалимо			
	Croatian	Vrlo zapaljivo			
	Czech	Vysoce hořlavé			
	Danish	Meget brandfarlig			
	Dutch	Licht ontvlambaar			
	Estonian	Kergesti süttiv			
	Finnish	Hyvin helposti syttyvä			
	French	Extrêmement inflammable			
	German	Entzündlich			
	Greek	Εξαιρετικά εύφλεκτο			
	Hungarian	Nagy mértékben gyúlékony			
	Icelandic	Mjög eldfimt			
	Italian	Altamente infiammabile			
	Latvian	Viegli uzliesmojošs			
	Lithuanian	Labai degus			
	Macedonian	Високо запаливо			
	Norwegian	Svært brannfarlig			
	Polish	Produkt wysoce łatwopalny			
	Portuguese	Altamente inflamável			
	Romanian	Extrem de inflamabil			
Serbian	Lako zapaljivo				
Slovak	Vysoko horľavý				
Slovenian	Visoko vnetljivo				
Spanish	Altamente inflamable				
Swedish	Mycket brandfarligt				
Turkish	Yüksek seviyede yanıcı				
	American English	Absorption capability	Indicates an absorption capability of medical device	N/A	N/A
	Bosnian	Sposobnost apsorpcije			
	British English	Absorption capability			
	Bulgarian	Абсорбционна способност			
	Croatian	Sposobnost apsorpcije			
	Czech	Absorpční kapacita			
	Danish	Indikation på absorberingsevne			
	Dutch	Absorberingsvermogen			
	Estonian	Imenduvus			
	Finnish	Imukyky			
	French	Capacité d'absorption			
	German	Absorptionsvermögen			
	Greek	Ικανότητα διαχείρισης του εξιδρώματος			
	Hungarian	Nedvszívó képesség			
	Icelandic	Ísogsgeta			
	Italian	Capacità di assorbimento			
	Latvian	Absorbcijas spēja			
	Lithuanian	Sugerties geba			
	Macedonian	Способност за апсорбирање			
	Norwegian	Absorpsjonsevne			
	Polish	Zdolność absorpcyjna			
	Portuguese	Capacidade de absorção			
	Romanian	Capacitate de absorbtie			
Serbian	Sposobnost apsorpcije				
Slovak	Kapacita absorpcie				
Slovenian	Zmožnost vpijanja				
Spanish	Capacidad de absorción				
Swedish	Absorptionsförmåga				
Turkish	Emilim kapasitesi				
	American English	Absorption capability	Indicates an absorption capability of medical device	N/A	N/A
	Bosnian	Sposobnost apsorpcije			
	British English	Absorption capability			
	Bulgarian	Абсорбционна способност			
	Croatian	Sposobnost apsorpcije			
	Czech	Absorpční kapacita			
	Danish	Indikation på absorberingsevne			
	Dutch	Absorberingsvermogen			
	Estonian	Imenduvus			
	Finnish	Imukyky			
	French	Capacité d'absorption			
	German	Absorptionsvermögen			
	Greek	Ικανότητα διαχείρισης του εξιδρώματος			
	Hungarian	Nedvszívó képesség			
	Icelandic	Ísogsgeta			
	Italian	Capacità di assorbimento			
	Latvian	Absorbcijas spēja			
	Lithuanian	Sugerties geba			
	Macedonian	Способност за апсорбирање			
	Norwegian	Absorpsjonsevne			
	Polish	Zdolność absorpcyjna			
	Portuguese	Capacidade de absorção			
	Romanian	Capacitate de absorbtie			
Serbian	Sposobnost apsorpcije				
Slovak	Kapacita absorpcie				
Slovenian	Zmožnost vpijanja				
Spanish	Capacidad de absorción				
Swedish	Absorptionsförmåga				
Turkish	Emilim kapasitesi				




	American English	Absorption capability	Indicates an absorption capability of medical device	N/A	N/A
	Bosnian	Sposobnost apsorpcije			
	British English	Absorption capability			
	Bulgarian	Абсорбционна способност			
	Croatian	Sposobnost apsorpcije			
	Czech	Absorpční kapacita			
	Danish	Indikation på absorberingsevne			
	Dutch	Absorberingsvermogen			
	Estonian	Imenduvus			
	Finnish	Imukyky			
	French	Capacité d'absorption			
	German	Absorptionsvermögen			
	Greek	Ικανότητα διαχείρισης του εξιδρώματος			
	Hungarian	Nedvszívó képesség			
	Icelandic	Ísogsgeta			
	Italian	Capacità di assorbimento			
	Latvian	Absorbcijas spēja			
	Lithuanian	Sugerties geba			
	Macedonian	Способност за апсорбирање			
	Norwegian	Absorpsjonsevne			
	Polish	Zdolność absorpcyjna			
	Portuguese	Capacidade de absorção			
	Romanian	Capacitate de absorbtie			
Serbian	Sposobnost apsorpcije				
Slovak	Kapacita absorpcie				
Slovenian	Zmožnost vpijanja				
Spanish	Capacidad de absorción				
Swedish	Absorptionsförmåga				
Turkish	Emilim kapasitesi				

Sterility




	American English	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.2.3.
	Bosnian	Sterilizovano etilen-oksidom			
	British English	Sterilized using ethylene oxide			
	Bulgarian	Стерилизирано с етиленов оксид			
	Croatian	Sterilizirano pomoću etilen oksida			
	Czech	Sterilizováno Ethylene Oxidem			
	Danish	Steriliseret ved brug af ethylenoxid			
	Dutch	Gesteriliseerd met Ethyleen Oxyde			
	Estonian	Steriliseeritud etüleenoksiidiga			
	Finnish	Steriloitu etyleenioksidilla			
	French	Stérilisé à l'oxyde d'éthylène			
	German	Mit Ethylenoxid sterilisiert			
	Greek	Αποστειρωμένο χρησιμοποιώντας οξείδιο του αιθυλενίου			
	Hungarian	Etilén-oxidál sterilizálva			
	Icelandic	Sæft með því að nota etýlenoxíð			
	Italian	Sterilizzato ad Ossido di Etilene			
	Latvian	Sterilizēts, izmantojot etilēnoksidu			
	Lithuanian	Sterilizuota etileno oksidu			
	Macedonian	Стерилизирано со етилен оксид			
	Norwegian	Sterilisert med Ethylene Oxide			
	Polish	Sterylizowany tlenkiem etylenu			
	Portuguese	Esterilizado usando óxido de etileno			
	Romanian	Sterilizat cu oxid de etilenă			
Serbian	Sterilisano etilen-oksidom				
Slovak	Sterilizované Ethylene Oxide				
Slovenian	Sterilizirano z etilenoksidom				
Spanish	Esterilización con óxido de etileno				
Swedish	Steriliserad med Ethylene Oxide				
Turkish	Etilen oksit kullanılarak sterilize edilmiştir				




	American English	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.2.4.
	Bosnian	Sterilizovano zračenjem			
	British English	Sterilized using irradiation			
	Bulgarian	Стерилизирано с радиация			
	Croatian	Sterilizirano zračenjem			
	Czech	Sterilizováno zářením			
	Danish	Steriliseret med bestråling			
	Dutch	Gesteriliseerd d.m.v. bestraling			
	Estonian	Steriliseeritud kiirgusega			
	Finnish	Steriloitu säteilyttämällä			
	French	Stérilisé par irradiation			
	German	Sterilisiert durch Bestrahlung			
	Greek	Αποστειρωμένο με τη χρήση ακτινοβολίας			
	Hungarian	Besugárzással sterilizálva			
	Icelandic	Sæft með geislun			
	Italian	Sterilizzato per irradiazione			
	Latvian	Sterilizēts, izmantojot starojumu			
	Lithuanian	Sterilizuota spinduliuote			
	Macedonian	Стерилизирано со употреба на зрачење			
	Norwegian	Sterilisert med stråling			
	Polish	Sterylizowany przez napromienianie			
	Portuguese	Esterilizado usando radiação			
	Romanian	Sterilizat prin iradiere			
Serbian	Sterilisano zračenjem				
Slovak	Sterilizované žiarením				
Slovenian	Izdelek je steriliziran z obsevanjem				
Spanish	Esterilizado con radiación				
Swedish	Steriliserad med strålning				
Turkish	İşinlama yöntemi ile steril edilmiştir				

	<table border="1"> <tbody> <tr><td>American English</td><td>Sterilized using steam or dry heat</td></tr> <tr><td>Bosnian</td><td>Sterilizovano parom ili suhim vrućim vazduhom</td></tr> <tr><td>British English</td><td>Sterilized using steam or dry heat</td></tr> <tr><td>Bulgarian</td><td>Стерилизирано с пара или суха топлина</td></tr> <tr><td>Croatian</td><td>Sterilizirano parom ili suhom toplinom</td></tr> <tr><td>Czech</td><td>Sterilizováno parou nebo suchým teplem</td></tr> <tr><td>Danish</td><td>Steriliseret vha. damp eller tør varme</td></tr> <tr><td>Dutch</td><td>Gesteriliseerd met stoom of droge hitte</td></tr> <tr><td>Estonian</td><td>Steriliseeritud auru või kuiva kuumusega</td></tr> <tr><td>Finnish</td><td>Steriloitu käyttäen höyryä tai kuumaa ilmaa</td></tr> <tr><td>French</td><td>Stérilisé à l'aide de vapeur ou de chaleur sèche</td></tr> <tr><td>German</td><td>Mit Dampf oder Trockenhitze sterilisiert</td></tr> <tr><td>Greek</td><td>Αποστειρωμένο με ατμό ή ξηρό αέρα</td></tr> <tr><td>Hungarian</td><td>Gőzzel vagy száraz hővel sterilizálva</td></tr> <tr><td>Icelandic</td><td>Sæft með því að nota gufu eða þurrán hita</td></tr> <tr><td>Italian</td><td>Sterilizzato a vapore o mediante calore secco</td></tr> <tr><td>Latvian</td><td>Sterilizēts, izmantojot tvaika vai sauso karšēšanu</td></tr> <tr><td>Lithuanian</td><td>Sterilizuota garu arba sausu karščiū</td></tr> <tr><td>Macedonian</td><td>Стерилизирано со пара или сува топлина</td></tr> <tr><td>Norwegian</td><td>Sterilisert med damp eller tørr varme</td></tr> <tr><td>Polish</td><td>Wyjałowione parą wodną lub gorącym suchym powietrzem</td></tr> <tr><td>Portuguese</td><td>Esterilizado através de vapor ou calor seco</td></tr> <tr><td>Romanian</td><td>Sterilizat cu aburi sau căldură uscată</td></tr> <tr><td>Serbian</td><td>Sterilisano parom ili suvom toplinom</td></tr> <tr><td>Slovak</td><td>Sterilizované pomocou pary alebo suchého tepla</td></tr> <tr><td>Slovenian</td><td>Sterilizirano s paro ali suho toploto</td></tr> <tr><td>Spanish</td><td>Esterilizado con vapor o calor seco</td></tr> <tr><td>Swedish</td><td>Steriliserad med ånga eller varmluft</td></tr> <tr><td>Turkish</td><td>Buhar veya kuru ısı kullanılarak sterilize edilmiştir</td></tr> </tbody> </table>	American English	Sterilized using steam or dry heat	Bosnian	Sterilizovano parom ili suhim vrućim vazduhom	British English	Sterilized using steam or dry heat	Bulgarian	Стерилизирано с пара или суха топлина	Croatian	Sterilizirano parom ili suhom toplinom	Czech	Sterilizováno parou nebo suchým teplem	Danish	Steriliseret vha. damp eller tør varme	Dutch	Gesteriliseerd met stoom of droge hitte	Estonian	Steriliseeritud auru või kuiva kuumusega	Finnish	Steriloitu käyttäen höyryä tai kuumaa ilmaa	French	Stérilisé à l'aide de vapeur ou de chaleur sèche	German	Mit Dampf oder Trockenhitze sterilisiert	Greek	Αποστειρωμένο με ατμό ή ξηρό αέρα	Hungarian	Gőzzel vagy száraz hővel sterilizálva	Icelandic	Sæft með því að nota gufu eða þurrán hita	Italian	Sterilizzato a vapore o mediante calore secco	Latvian	Sterilizēts, izmantojot tvaika vai sauso karšēšanu	Lithuanian	Sterilizuota garu arba sausu karščiū	Macedonian	Стерилизирано со пара или сува топлина	Norwegian	Sterilisert med damp eller tørr varme	Polish	Wyjałowione parą wodną lub gorącym suchym powietrzem	Portuguese	Esterilizado através de vapor ou calor seco	Romanian	Sterilizat cu aburi sau căldură uscată	Serbian	Sterilisano parom ili suvom toplinom	Slovak	Sterilizované pomocou pary alebo suchého tepla	Slovenian	Sterilizirano s paro ali suho toploto	Spanish	Esterilizado con vapor o calor seco	Swedish	Steriliserad med ånga eller varmluft	Turkish	Buhar veya kuru ısı kullanılarak sterilize edilmiştir	<p>Indicates a medical device that has been sterilized using steam or dry heat</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.2.5.</p>
American English	Sterilized using steam or dry heat																																																													
Bosnian	Sterilizovano parom ili suhim vrućim vazduhom																																																													
British English	Sterilized using steam or dry heat																																																													
Bulgarian	Стерилизирано с пара или суха топлина																																																													
Croatian	Sterilizirano parom ili suhom toplinom																																																													
Czech	Sterilizováno parou nebo suchým teplem																																																													
Danish	Steriliseret vha. damp eller tør varme																																																													
Dutch	Gesteriliseerd met stoom of droge hitte																																																													
Estonian	Steriliseeritud auru või kuiva kuumusega																																																													
Finnish	Steriloitu käyttäen höyryä tai kuumaa ilmaa																																																													
French	Stérilisé à l'aide de vapeur ou de chaleur sèche																																																													
German	Mit Dampf oder Trockenhitze sterilisiert																																																													
Greek	Αποστειρωμένο με ατμό ή ξηρό αέρα																																																													
Hungarian	Gőzzel vagy száraz hővel sterilizálva																																																													
Icelandic	Sæft með því að nota gufu eða þurrán hita																																																													
Italian	Sterilizzato a vapore o mediante calore secco																																																													
Latvian	Sterilizēts, izmantojot tvaika vai sauso karšēšanu																																																													
Lithuanian	Sterilizuota garu arba sausu karščiū																																																													
Macedonian	Стерилизирано со пара или сува топлина																																																													
Norwegian	Sterilisert med damp eller tørr varme																																																													
Polish	Wyjałowione parą wodną lub gorącym suchym powietrzem																																																													
Portuguese	Esterilizado através de vapor ou calor seco																																																													
Romanian	Sterilizat cu aburi sau căldură uscată																																																													
Serbian	Sterilisano parom ili suvom toplinom																																																													
Slovak	Sterilizované pomocou pary alebo suchého tepla																																																													
Slovenian	Sterilizirano s paro ali suho toploto																																																													
Spanish	Esterilizado con vapor o calor seco																																																													
Swedish	Steriliserad med ånga eller varmluft																																																													
Turkish	Buhar veya kuru ısı kullanılarak sterilize edilmiştir																																																													
	<table border="1"> <tbody> <tr><td>American English</td><td>Non-sterile</td></tr> <tr><td>Bosnian</td><td>Nije sterilno</td></tr> <tr><td>British English</td><td>Non-sterile</td></tr> <tr><td>Bulgarian</td><td>Нестерилно</td></tr> <tr><td>Croatian</td><td>Nesterilno</td></tr> <tr><td>Czech</td><td>Nesterilní</td></tr> <tr><td>Danish</td><td>Usteril</td></tr> <tr><td>Dutch</td><td>Niet-steriel</td></tr> <tr><td>Estonian</td><td>Mittesteriilne</td></tr> <tr><td>Finnish</td><td>Steriloimaton</td></tr> <tr><td>French</td><td>Non stérile</td></tr> <tr><td>German</td><td>Unsteril</td></tr> <tr><td>Greek</td><td>Μη αποστειρωμένο</td></tr> <tr><td>Hungarian</td><td>Nem steril</td></tr> <tr><td>Icelandic</td><td>EKKI dauðhreinsað</td></tr> <tr><td>Italian</td><td>Non sterile</td></tr> <tr><td>Latvian</td><td>Nav sterilis</td></tr> <tr><td>Lithuanian</td><td>Nesterilu</td></tr> <tr><td>Macedonian</td><td>Не е стерилно</td></tr> <tr><td>Norwegian</td><td>Usteril</td></tr> <tr><td>Polish</td><td>Niesterylony</td></tr> <tr><td>Portuguese</td><td>Não estéril</td></tr> <tr><td>Romanian</td><td>Nesteril</td></tr> <tr><td>Serbian</td><td>Nesterilno</td></tr> <tr><td>Slovak</td><td>Nesterilné</td></tr> <tr><td>Slovenian</td><td>Nesterilno</td></tr> <tr><td>Spanish</td><td>No estéril</td></tr> <tr><td>Swedish</td><td>Icke-steril</td></tr> <tr><td>Turkish</td><td>Steril değildir</td></tr> </tbody> </table>	American English	Non-sterile	Bosnian	Nije sterilno	British English	Non-sterile	Bulgarian	Нестерилно	Croatian	Nesterilno	Czech	Nesterilní	Danish	Usteril	Dutch	Niet-steriel	Estonian	Mittesteriilne	Finnish	Steriloimaton	French	Non stérile	German	Unsteril	Greek	Μη αποστειρωμένο	Hungarian	Nem steril	Icelandic	EKKI dauðhreinsað	Italian	Non sterile	Latvian	Nav sterilis	Lithuanian	Nesterilu	Macedonian	Не е стерилно	Norwegian	Usteril	Polish	Niesterylony	Portuguese	Não estéril	Romanian	Nesteril	Serbian	Nesterilno	Slovak	Nesterilné	Slovenian	Nesterilno	Spanish	No estéril	Swedish	Icke-steril	Turkish	Steril değildir	<p>Indicates a medical device that has not been subjected to a sterilization process</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.2.7.</p>
American English	Non-sterile																																																													
Bosnian	Nije sterilno																																																													
British English	Non-sterile																																																													
Bulgarian	Нестерилно																																																													
Croatian	Nesterilno																																																													
Czech	Nesterilní																																																													
Danish	Usteril																																																													
Dutch	Niet-steriel																																																													
Estonian	Mittesteriilne																																																													
Finnish	Steriloimaton																																																													
French	Non stérile																																																													
German	Unsteril																																																													
Greek	Μη αποστειρωμένο																																																													
Hungarian	Nem steril																																																													
Icelandic	EKKI dauðhreinsað																																																													
Italian	Non sterile																																																													
Latvian	Nav sterilis																																																													
Lithuanian	Nesterilu																																																													
Macedonian	Не е стерилно																																																													
Norwegian	Usteril																																																													
Polish	Niesterylony																																																													
Portuguese	Não estéril																																																													
Romanian	Nesteril																																																													
Serbian	Nesterilno																																																													
Slovak	Nesterilné																																																													
Slovenian	Nesterilno																																																													
Spanish	No estéril																																																													
Swedish	Icke-steril																																																													
Turkish	Steril değildir																																																													
	<table border="1"> <tbody> <tr><td>American English</td><td>Do not use if package is damaged and consult instructions for use</td></tr> <tr><td>Bosnian</td><td>Nemojte koristiti ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu</td></tr> <tr><td>British English</td><td>Do not use if package is damaged and consult instructions for use</td></tr> <tr><td>Bulgarian</td><td>Да не се използва, ако опаковката е повредена и да се направи справка с инструкциите за употреба</td></tr> <tr><td>Croatian</td><td>Nemojte koristiti ako je pakiranje oštećeno te proučite upute za uporabu</td></tr> <tr><td>Czech</td><td>Pokud je obal poškozený, nepoužívejte a přečtěte si návod k použití</td></tr> <tr><td>Danish</td><td>Ma ikke bruges, hvis emballagen er beskadiget, og se brugsanvisningen</td></tr> <tr><td>Dutch</td><td>Niet gebruiken als de verpakking beschadigd is en raadpleeg de gebruiksaanwijzing</td></tr> <tr><td>Estonian</td><td>Mitte kasutada, kui pakend on kahjustatud, ja lugeda kasutusjuhendit</td></tr> <tr><td>Finnish</td><td>Älä käytä jos pakkaus on vaurioitunut ja lue käyttöohjeet ennen käyttöä</td></tr> <tr><td>French</td><td>Ne pas utiliser si l'emballage est endommagé et consulter le mode d'emploi</td></tr> <tr><td>German</td><td>Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten</td></tr> <tr><td>Greek</td><td>Να μη χρησιμοποιείται αν η συσκευασία έχει υποστεί ζημιά και ανατρέξτε στις οδηγίες χρήσης</td></tr> <tr><td>Hungarian</td><td>Ne használja, ha a csomag sérült, és olvassa el a használati utasítást</td></tr> <tr><td>Icelandic</td><td>Skal ekki nota ef innsigli þakningar er rofið og skoðið leiðbeiningar fyrir notkun</td></tr> <tr><td>Italian</td><td>Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso</td></tr> <tr><td>Latvian</td><td>Nelietot, ja iepakojums bojāts, un skatīt lietošanas pamācību</td></tr> <tr><td>Lithuanian</td><td>Nenaudoti, jei pakuotė pažeista, ir skaityti naudojimo instrukciją</td></tr> <tr><td>Macedonian</td><td>Да не се користи ако пакувањето е оштетено и да се проверат упатствата за употреба</td></tr> <tr><td>Norwegian</td><td>Ikke bruk hvis pakningen er skadd, og se bruksanvisningen</td></tr> <tr><td>Polish</td><td>Nie używać, jeżeli opakowanie jest uszkodzone i zapoznaj się z instrukcją używania</td></tr> <tr><td>Portuguese</td><td>Não utilize se a embalagem estiver danificada e consulte as instruções de utilização</td></tr> <tr><td>Romanian</td><td>A nu se utiliza dacă ambalajul este deteriorat și a se consulta instrucțiunile de utilizare</td></tr> <tr><td>Serbian</td><td>Nemojte koristiti ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu</td></tr> <tr><td>Slovak</td><td>Nepoužívajte, ak je obal poškodený, alebo si prečítajte pokyny na používanie</td></tr> <tr><td>Slovenian</td><td>Ne uporabljajte, če je ovojnina poškodovana, in upoštevajte navodila za uporabo</td></tr> <tr><td>Spanish</td><td>No utilizar si el envase está dañado y consultar las instrucciones de uso</td></tr> <tr><td>Swedish</td><td>Använd inte om förpackningen är skadad och läs bruksanvisningen</td></tr> <tr><td>Turkish</td><td>Ambalaj hasarlıysa kullanmayın ve kullanim talimatina başvurun</td></tr> </tbody> </table>	American English	Do not use if package is damaged and consult instructions for use	Bosnian	Nemojte koristiti ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	British English	Do not use if package is damaged and consult instructions for use	Bulgarian	Да не се използва, ако опаковката е повредена и да се направи справка с инструкциите за употреба	Croatian	Nemojte koristiti ako je pakiranje oštećeno te proučite upute za uporabu	Czech	Pokud je obal poškozený, nepoužívejte a přečtěte si návod k použití	Danish	Ma ikke bruges, hvis emballagen er beskadiget, og se brugsanvisningen	Dutch	Niet gebruiken als de verpakking beschadigd is en raadpleeg de gebruiksaanwijzing	Estonian	Mitte kasutada, kui pakend on kahjustatud, ja lugeda kasutusjuhendit	Finnish	Älä käytä jos pakkaus on vaurioitunut ja lue käyttöohjeet ennen käyttöä	French	Ne pas utiliser si l'emballage est endommagé et consulter le mode d'emploi	German	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	Greek	Να μη χρησιμοποιείται αν η συσκευασία έχει υποστεί ζημιά και ανατρέξτε στις οδηγίες χρήσης	Hungarian	Ne használja, ha a csomag sérült, és olvassa el a használati utasítást	Icelandic	Skal ekki nota ef innsigli þakningar er rofið og skoðið leiðbeiningar fyrir notkun	Italian	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso	Latvian	Nelietot, ja iepakojums bojāts, un skatīt lietošanas pamācību	Lithuanian	Nenaudoti, jei pakuotė pažeista, ir skaityti naudojimo instrukciją	Macedonian	Да не се користи ако пакувањето е оштетено и да се проверат упатствата за употреба	Norwegian	Ikke bruk hvis pakningen er skadd, og se bruksanvisningen	Polish	Nie używać, jeżeli opakowanie jest uszkodzone i zapoznaj się z instrukcją używania	Portuguese	Não utilize se a embalagem estiver danificada e consulte as instruções de utilização	Romanian	A nu se utiliza dacă ambalajul este deteriorat și a se consulta instrucțiunile de utilizare	Serbian	Nemojte koristiti ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Slovak	Nepoužívajte, ak je obal poškodený, alebo si prečítajte pokyny na používanie	Slovenian	Ne uporabljajte, če je ovojnina poškodovana, in upoštevajte navodila za uporabo	Spanish	No utilizar si el envase está dañado y consultar las instrucciones de uso	Swedish	Använd inte om förpackningen är skadad och läs bruksanvisningen	Turkish	Ambalaj hasarlıysa kullanmayın ve kullanim talimatina başvurun	<p>Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.2.8.</p>
American English	Do not use if package is damaged and consult instructions for use																																																													
Bosnian	Nemojte koristiti ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu																																																													
British English	Do not use if package is damaged and consult instructions for use																																																													
Bulgarian	Да не се използва, ако опаковката е повредена и да се направи справка с инструкциите за употреба																																																													
Croatian	Nemojte koristiti ako je pakiranje oštećeno te proučite upute za uporabu																																																													
Czech	Pokud je obal poškozený, nepoužívejte a přečtěte si návod k použití																																																													
Danish	Ma ikke bruges, hvis emballagen er beskadiget, og se brugsanvisningen																																																													
Dutch	Niet gebruiken als de verpakking beschadigd is en raadpleeg de gebruiksaanwijzing																																																													
Estonian	Mitte kasutada, kui pakend on kahjustatud, ja lugeda kasutusjuhendit																																																													
Finnish	Älä käytä jos pakkaus on vaurioitunut ja lue käyttöohjeet ennen käyttöä																																																													
French	Ne pas utiliser si l'emballage est endommagé et consulter le mode d'emploi																																																													
German	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten																																																													
Greek	Να μη χρησιμοποιείται αν η συσκευασία έχει υποστεί ζημιά και ανατρέξτε στις οδηγίες χρήσης																																																													
Hungarian	Ne használja, ha a csomag sérült, és olvassa el a használati utasítást																																																													
Icelandic	Skal ekki nota ef innsigli þakningar er rofið og skoðið leiðbeiningar fyrir notkun																																																													
Italian	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso																																																													
Latvian	Nelietot, ja iepakojums bojāts, un skatīt lietošanas pamācību																																																													
Lithuanian	Nenaudoti, jei pakuotė pažeista, ir skaityti naudojimo instrukciją																																																													
Macedonian	Да не се користи ако пакувањето е оштетено и да се проверат упатствата за употреба																																																													
Norwegian	Ikke bruk hvis pakningen er skadd, og se bruksanvisningen																																																													
Polish	Nie używać, jeżeli opakowanie jest uszkodzone i zapoznaj się z instrukcją używania																																																													
Portuguese	Não utilize se a embalagem estiver danificada e consulte as instruções de utilização																																																													
Romanian	A nu se utiliza dacă ambalajul este deteriorat și a se consulta instrucțiunile de utilizare																																																													
Serbian	Nemojte koristiti ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu																																																													
Slovak	Nepoužívajte, ak je obal poškodený, alebo si prečítajte pokyny na používanie																																																													
Slovenian	Ne uporabljajte, če je ovojnina poškodovana, in upoštevajte navodila za uporabo																																																													
Spanish	No utilizar si el envase está dañado y consultar las instrucciones de uso																																																													
Swedish	Använd inte om förpackningen är skadad och läs bruksanvisningen																																																													
Turkish	Ambalaj hasarlıysa kullanmayın ve kullanim talimatina başvurun																																																													

	American English	Single sterile barrier system	Indicates a single sterile barrier system	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.2.11.
	Bosnian	Sistem jednostruke sterilne barijere			
	British English	Single sterile barrier system			
	Bulgarian	Единична стерилна бариерна система			
	Croatian	Sustav jednostruke sterilne barijere			
	Czech	Jednodilný sterilní bariérový systém			
	Danish	Indikerer et sterilt beskyttelsessystem til engangsbrug			
	Dutch	Enkel steriel barrièresysteem			
	Estonian	Ühekordne steriilne kaitsepakend			
	Finnish	Yksittäinen steriili suojajärjestelmä			
	French	Indique un système à barrière unique stérile			
	German	System mit nur einer Sterilbarriere			
	Greek	Μοναδικό αποστειρωμένο σύστημα συγκράτησης			
	Hungarian	Egyszeres steril védőzárendszer			
	Icelandic	Stakt saft hindrandi kerfi			
	Italian	Indica un sistema di barriera sterile singola			
	Latvian	Viena sterila barjeras sistēma			
	Lithuanian	Viengubo sterilaus barjero sistema			
	Macedonian	Систем со единична стерилна бариера			
	Norwegian	System med én steril barriere			
	Polish	System pojedynczej bariery jalowej			
Portuguese	Sistema de barreira estéril única				
Romanian	Sistem cu barieră sterilă unică				
Serbian	Sistem pojedinačne sterilne barijere				
Slovak	Systém jednej sterilnej bariéry				
Slovenian	Sistem enojne sterilne pregrade				
Spanish	Indica un único sistema de barrera estéril				
Swedish	Enkelt, sterilit barriärsystem				
Turkish	Tek steril bariyer sistemi				
	American English	Single sterile barrier system	Indicates a single sterile barrier system	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.2.11.
	Bosnian	Sistem jednostruke sterilne barijere			
	British English	Single sterile barrier system			
	Bulgarian	Единична стерилна бариерна система			
	Croatian	Sustav jednostruke sterilne barijere			
	Czech	Jednodilný sterilní bariérový systém			
	Danish	Indikerer et sterilt beskyttelsessystem til engangsbrug			
	Dutch	Enkel steriel barrièresysteem			
	Estonian	Ühekordne steriilne kaitsepakend			
	Finnish	Yksittäinen steriili suojajärjestelmä			
	French	Indique un système à barrière unique stérile			
	German	System mit nur einer Sterilbarriere			
	Greek	Μοναδικό αποστειρωμένο σύστημα συγκράτησης			
	Hungarian	Egyszeres steril védőzárendszer			
	Icelandic	Stakt saft hindrandi kerfi			
	Italian	Indica un sistema di barriera sterile singola			
	Latvian	Viena sterila barjeras sistēma			
	Lithuanian	Viengubo sterilaus barjero sistema			
	Macedonian	Систем со единична стерилна бариера			
	Norwegian	System med én steril barriere			
	Polish	System pojedynczej bariery jalowej			
Portuguese	Sistema de barreira estéril única				
Romanian	Sistem cu barieră sterilă unică				
Serbian	Sistem pojedinačne sterilne barijere				
Slovak	Systém jednej sterilnej bariéry				
Slovenian	Sistem enojne sterilne pregrade				
Spanish	Indica un único sistema de barrera estéril				
Swedish	Enkelt, sterilit barriärsystem				
Turkish	Tek steril bariyer sistemi				
	American English	Single sterile barrier system	Indicates a single sterile barrier system	ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	5.2.11.
	Bosnian	Sistem jednostruke sterilne barijere			
	British English	Single sterile barrier system			
	Bulgarian	Единична стерилна бариерна система			
	Croatian	Sustav jednostruke sterilne barijere			
	Czech	Jednodilný sterilní bariérový systém			
	Danish	Indikerer et sterilt beskyttelsessystem til engangsbrug			
	Dutch	Enkel steriel barrièresysteem			
	Estonian	Ühekordne steriilne kaitsepakend			
	Finnish	Yksittäinen steriili suojajärjestelmä			
	French	Indique un système à barrière unique stérile			
	German	System mit nur einer Sterilbarriere			
	Greek	Μοναδικό αποστειρωμένο σύστημα συγκράτησης			
	Hungarian	Egyszeres steril védőzárendszer			
	Icelandic	Stakt saft hindrandi kerfi			
	Italian	Indica un sistema di barriera sterile singola			
	Latvian	Viena sterila barjeras sistēma			
	Lithuanian	Viengubo sterilaus barjero sistema			
	Macedonian	Систем со единична стерилна бариера			
	Norwegian	System med én steril barriere			
	Polish	System pojedynczej bariery jalowej			
Portuguese	Sistema de barreira estéril única				
Romanian	Sistem cu barieră sterilă unică				
Serbian	Sistem pojedinačne sterilne barijere				
Slovak	Systém jednej sterilnej bariéry				
Slovenian	Sistem enojne sterilne pregrade				
Spanish	Indica un único sistema de barrera estéril				
Swedish	Enkelt, sterilit barriärsystem				
Turkish	Tek steril bariyer sistemi				

Storage

	<p>American English Keep away from sunlight</p> <p>Bosnian Držite dalje od sunčeve svjetlosti</p> <p>British English Keep away from sunlight</p> <p>Bulgarian Да се пазит от слънчева светлина</p> <p>Croatian Držati podalje od sunčeve svjetlosti</p> <p>Czech Uchovávat mimo sluneční záření</p> <p>Danish Bør ikke udsættes for sollys</p> <p>Dutch Uit de buurt van zonlicht houden</p> <p>Estonian Hoida päkesevalgusest eemal.</p> <p>Finnish Suojattava suoralta auringonvalolta</p> <p>French Conserver à l'abri de la lumière</p> <p>German Vor Sonnenlicht schützen</p> <p>Greek Αποφύγετε την απευθείας έκθεση στον ήλιο</p> <p>Hungarian Napfénytől elzárva tartandó</p> <p>Icelandic Geymdu fjarri sólarljósi</p> <p>Italian Conservare lontano dai raggi solari</p> <p>Latvian Sargāt no saules gaismas</p> <p>Lithuanian Saugoti nuo saulės šviesos</p> <p>Macedonian Чувајте го подалеку од сончева светлина</p> <p>Norwegian Holdes unna sollys</p> <p>Polish Przechowywać z dala od światła słonecznego</p> <p>Portuguese Manter afastado da luz solar</p> <p>Romanian A se feri de lumina soarelui</p> <p>Serbian Držati podalje od sunčeve svetlosti</p> <p>Slovak Chráňte pred slnečným svetlom</p> <p>Slovenian Izdelka ne izpostavljajte soncu</p> <p>Spanish Mantener alejado de la luz solar</p> <p>Swedish Skyddas från solljus</p> <p>Turkish Güneş ışığından uzak tutun</p>	<p>Indicates a medical device that needs protection from light sources</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.3.2.</p>
	<p>American English Lower limit of temperature</p> <p>Bosnian Donja temperaturna granica</p> <p>British English Lower limit of temperature</p> <p>Bulgarian Долна граница на температурата</p> <p>Croatian Donje ograničenje temperature</p> <p>Czech Dolní teplotní limit</p> <p>Danish Nedre temperaturgrænse</p> <p>Dutch Ondergrens temperatuur</p> <p>Estonian Temperatuuri alampiir</p> <p>Finnish Lämpötilan alaraja</p> <p>French Limite inférieure de température</p> <p>German Unterer Temperaturgrenzwert</p> <p>Greek Κατώτατο όριο θερμοκρασίας</p> <p>Hungarian Alsó hőmérsékleti határ</p> <p>Icelandic Lægi mörk hitastigs</p> <p>Italian Limite inferiore di temperatura</p> <p>Latvian Temperatūras zemākā robeža</p> <p>Lithuanian Apatinė temperatūros riba</p> <p>Macedonian Долна температурна граница</p> <p>Norwegian Nedre temperaturgrænse</p> <p>Polish Dolny limit temperatury</p> <p>Portuguese Limite inferior de temperatura</p> <p>Romanian Limită inferioară de temperatură</p> <p>Serbian Donja granica temperature</p> <p>Slovak Najnižšia prípustná teplota</p> <p>Slovenian Spodnja omejitve temperature</p> <p>Spanish Limite inferior de temperatura</p> <p>Swedish Nedre temperaturgräns</p> <p>Turkish Sıcaklık alt sınırı</p>	<p>Indicates the lower limit of temperature to which the medical device can be safely exposed</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.3.5.</p>
	<p>American English Upper limit of temperature</p> <p>Bosnian Gornja temperaturna granica</p> <p>British English Upper limit of temperature</p> <p>Bulgarian Горна граница на температурата</p> <p>Croatian Gornje ograničenje temperature</p> <p>Czech Horní teplotní limit</p> <p>Danish Øvre temperaturgrænse</p> <p>Dutch Bovengrens temperatuur</p> <p>Estonian Temperatuuri ülempiir</p> <p>Finnish Lämpötilan yläraja</p> <p>French Limite supérieure de température</p> <p>German Oberer Temperaturgrenzwert</p> <p>Greek Ανώτατο όριο θερμοκρασίας</p> <p>Hungarian Felső hőmérsékleti határ</p> <p>Icelandic Efri mörk hitastigs</p> <p>Italian Limite superiore di temperatura</p> <p>Latvian Temperatūras augstākā robeža</p> <p>Lithuanian Viršutinė temperatūros riba</p> <p>Macedonian Горна температурна граница</p> <p>Norwegian Øvre temperaturgrænse</p> <p>Polish Górny limit temperatury</p> <p>Portuguese Limite superior de temperatura</p> <p>Romanian Limită superioară de temperatură</p> <p>Serbian Gornja granica temperature</p> <p>Slovak Najvyššia prípustná teplota</p> <p>Slovenian Zgornja omejitve temperature</p> <p>Spanish Limite superior de temperatura</p> <p>Swedish Övre temperaturgräns</p> <p>Turkish Sıcaklık üst sınırı</p>	<p>Indicates the upper limit of temperature to which the medical device can be safely exposed</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.3.6.</p>

	<table border="1"> <tr><td>American English</td><td>Temperature limit</td></tr> <tr><td>Bosnian</td><td>Temperaturna granica</td></tr> <tr><td>British English</td><td>Temperature limit</td></tr> <tr><td>Bulgarian</td><td>Температурна граница</td></tr> <tr><td>Croatian</td><td>Ograničenje temperature</td></tr> <tr><td>Czech</td><td>Teplotní limit</td></tr> <tr><td>Danish</td><td>Temperaturgrænse</td></tr> <tr><td>Dutch</td><td>Temperatuurgrens</td></tr> <tr><td>Estonian</td><td>Temperatuuripiir</td></tr> <tr><td>Finnish</td><td>Lämpötilaraja</td></tr> <tr><td>French</td><td>Limite de température</td></tr> <tr><td>German</td><td>Temperaturbegrenzung</td></tr> <tr><td>Greek</td><td>Όριο θερμοκρασίας</td></tr> <tr><td>Hungarian</td><td>Hőmérsékleti határ</td></tr> <tr><td>Icelandic</td><td>Hítastigsmörk</td></tr> <tr><td>Italian</td><td>Limite di temperatura</td></tr> <tr><td>Latvian</td><td>Temperatūras ierobežojums</td></tr> <tr><td>Lithuanian</td><td>Temperatūros ribos</td></tr> <tr><td>Macedonian</td><td>Температурна граница</td></tr> <tr><td>Norwegian</td><td>Temperaturgrense</td></tr> <tr><td>Polish</td><td>Limit temperatury</td></tr> <tr><td>Portuguese</td><td>Limite de temperatura</td></tr> <tr><td>Romanian</td><td>Limită de temperatură</td></tr> <tr><td>Serbian</td><td>Maksimalna temperatura</td></tr> <tr><td>Slovak</td><td>Teplotné obmedzenie</td></tr> <tr><td>Slovenian</td><td>Temperaturna omejitev</td></tr> <tr><td>Spanish</td><td>Limite de temperatura</td></tr> <tr><td>Swedish</td><td>Temperaturgräns</td></tr> <tr><td>Turkish</td><td>Sıcaklık sınırı</td></tr> </table>	American English	Temperature limit	Bosnian	Temperaturna granica	British English	Temperature limit	Bulgarian	Температурна граница	Croatian	Ograničenje temperature	Czech	Teplotní limit	Danish	Temperaturgrænse	Dutch	Temperatuurgrens	Estonian	Temperatuuripiir	Finnish	Lämpötilaraja	French	Limite de température	German	Temperaturbegrenzung	Greek	Όριο θερμοκρασίας	Hungarian	Hőmérsékleti határ	Icelandic	Hítastigsmörk	Italian	Limite di temperatura	Latvian	Temperatūras ierobežojums	Lithuanian	Temperatūros ribos	Macedonian	Температурна граница	Norwegian	Temperaturgrense	Polish	Limit temperatury	Portuguese	Limite de temperatura	Romanian	Limită de temperatură	Serbian	Maksimalna temperatura	Slovak	Teplotné obmedzenie	Slovenian	Temperaturna omejitev	Spanish	Limite de temperatura	Swedish	Temperaturgräns	Turkish	Sıcaklık sınırı	<p>Indicates the temperature limits to which the medical device can be safely exposed</p>	<p>ISO 15223-1:2020 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements</p>	<p>5.3.7.</p>
American English	Temperature limit																																																													
Bosnian	Temperaturna granica																																																													
British English	Temperature limit																																																													
Bulgarian	Температурна граница																																																													
Croatian	Ograničenje temperature																																																													
Czech	Teplotní limit																																																													
Danish	Temperaturgrænse																																																													
Dutch	Temperatuurgrens																																																													
Estonian	Temperatuuripiir																																																													
Finnish	Lämpötilaraja																																																													
French	Limite de température																																																													
German	Temperaturbegrenzung																																																													
Greek	Όριο θερμοκρασίας																																																													
Hungarian	Hőmérsékleti határ																																																													
Icelandic	Hítastigsmörk																																																													
Italian	Limite di temperatura																																																													
Latvian	Temperatūras ierobežojums																																																													
Lithuanian	Temperatūros ribos																																																													
Macedonian	Температурна граница																																																													
Norwegian	Temperaturgrense																																																													
Polish	Limit temperatury																																																													
Portuguese	Limite de temperatura																																																													
Romanian	Limită de temperatură																																																													
Serbian	Maksimalna temperatura																																																													
Slovak	Teplotné obmedzenie																																																													
Slovenian	Temperaturna omejitev																																																													
Spanish	Limite de temperatura																																																													
Swedish	Temperaturgräns																																																													
Turkish	Sıcaklık sınırı																																																													
	<table border="1"> <tr><td>American English</td><td>This way up</td></tr> <tr><td>Bosnian</td><td>Ova strana gore</td></tr> <tr><td>British English</td><td>This way up</td></tr> <tr><td>Bulgarian</td><td>Тази страна нагоре</td></tr> <tr><td>Croatian</td><td>Ovaj dio okrenuti prema gore</td></tr> <tr><td>Czech</td><td>Neklopit</td></tr> <tr><td>Danish</td><td>Denne side opad</td></tr> <tr><td>Dutch</td><td>Deze kant boven</td></tr> <tr><td>Estonian</td><td>See pool üles</td></tr> <tr><td>Finnish</td><td>Tämä puoli ylöspäin</td></tr> <tr><td>French</td><td>Haut</td></tr> <tr><td>German</td><td>Hier oben</td></tr> <tr><td>Greek</td><td>Προς τα πάνω</td></tr> <tr><td>Hungarian</td><td>Felfelé</td></tr> <tr><td>Icelandic</td><td>Þessi hlið upp</td></tr> <tr><td>Italian</td><td>Conservare in posizione verticale con le frecce rivolte verso l'alto</td></tr> <tr><td>Latvian</td><td>Šādi uz augšu</td></tr> <tr><td>Lithuanian</td><td>Šia puse į viršų</td></tr> <tr><td>Macedonian</td><td>Оваа страна нагоре</td></tr> <tr><td>Norwegian</td><td>Denne veien opp</td></tr> <tr><td>Polish</td><td>Tą stroną do góry</td></tr> <tr><td>Portuguese</td><td>Este lado para cima</td></tr> <tr><td>Romanian</td><td>Cu această parte în sus</td></tr> <tr><td>Serbian</td><td>Ovaj kraj nagore</td></tr> <tr><td>Slovak</td><td>Touto stranou nahor</td></tr> <tr><td>Slovenian</td><td>Ta stran navzgor</td></tr> <tr><td>Spanish</td><td>Este lado hacia arriba</td></tr> <tr><td>Swedish</td><td>Denna sida upp</td></tr> <tr><td>Turkish</td><td>Bu taraf yukarı</td></tr> </table>	American English	This way up	Bosnian	Ova strana gore	British English	This way up	Bulgarian	Тази страна нагоре	Croatian	Ovaj dio okrenuti prema gore	Czech	Neklopit	Danish	Denne side opad	Dutch	Deze kant boven	Estonian	See pool üles	Finnish	Tämä puoli ylöspäin	French	Haut	German	Hier oben	Greek	Προς τα πάνω	Hungarian	Felfelé	Icelandic	Þessi hlið upp	Italian	Conservare in posizione verticale con le frecce rivolte verso l'alto	Latvian	Šādi uz augšu	Lithuanian	Šia puse į viršų	Macedonian	Оваа страна нагоре	Norwegian	Denne veien opp	Polish	Tą stroną do góry	Portuguese	Este lado para cima	Romanian	Cu această parte în sus	Serbian	Ovaj kraj nagore	Slovak	Touto stranou nahor	Slovenian	Ta stran navzgor	Spanish	Este lado hacia arriba	Swedish	Denna sida upp	Turkish	Bu taraf yukarı	<p>To indicate correct upright position of the transport package</p>	<p>ISO 7000 Graphical symbols for use on equipment — Registered symbols</p>	<p>0623</p>
American English	This way up																																																													
Bosnian	Ova strana gore																																																													
British English	This way up																																																													
Bulgarian	Тази страна нагоре																																																													
Croatian	Ovaj dio okrenuti prema gore																																																													
Czech	Neklopit																																																													
Danish	Denne side opad																																																													
Dutch	Deze kant boven																																																													
Estonian	See pool üles																																																													
Finnish	Tämä puoli ylöspäin																																																													
French	Haut																																																													
German	Hier oben																																																													
Greek	Προς τα πάνω																																																													
Hungarian	Felfelé																																																													
Icelandic	Þessi hlið upp																																																													
Italian	Conservare in posizione verticale con le frecce rivolte verso l'alto																																																													
Latvian	Šādi uz augšu																																																													
Lithuanian	Šia puse į viršų																																																													
Macedonian	Оваа страна нагоре																																																													
Norwegian	Denne veien opp																																																													
Polish	Tą stroną do góry																																																													
Portuguese	Este lado para cima																																																													
Romanian	Cu această parte în sus																																																													
Serbian	Ovaj kraj nagore																																																													
Slovak	Touto stranou nahor																																																													
Slovenian	Ta stran navzgor																																																													
Spanish	Este lado hacia arriba																																																													
Swedish	Denna sida upp																																																													
Turkish	Bu taraf yukarı																																																													
	<table border="1"> <tr><td>American English</td><td>Recyclable packaging</td></tr> <tr><td>Bosnian</td><td>Pakovanje se može reciklirati</td></tr> <tr><td>British English</td><td>Recyclable packaging</td></tr> <tr><td>Bulgarian</td><td>Рециклируема опаковка</td></tr> <tr><td>Croatian</td><td>Pakiranje koje se može reciklirati</td></tr> <tr><td>Czech</td><td>Recyklovatelný obal</td></tr> <tr><td>Danish</td><td>Genanvendelig emballage</td></tr> <tr><td>Dutch</td><td>Recyclebare verpakking</td></tr> <tr><td>Estonian</td><td>Taaskasutatav pakend</td></tr> <tr><td>Finnish</td><td>Kierrätettävä pakkausmateriaali</td></tr> <tr><td>French</td><td>Emballage recyclable</td></tr> <tr><td>German</td><td>Recyclingfähige Verpackung</td></tr> <tr><td>Greek</td><td>Ανακυκλώσιμη συσκευασία</td></tr> <tr><td>Hungarian</td><td>Újrahasznosítható csomagolás</td></tr> <tr><td>Icelandic</td><td>Endurvinnanlegur umbúðir</td></tr> <tr><td>Italian</td><td>Confezione riciclabile</td></tr> <tr><td>Latvian</td><td>Pārstrādājams iepakojums</td></tr> <tr><td>Lithuanian</td><td>Perdirbama pakuotė</td></tr> <tr><td>Macedonian</td><td>Амбалажата може да се рециклира</td></tr> <tr><td>Norwegian</td><td>Resirkulerbar pakning</td></tr> <tr><td>Polish</td><td>Opakowanie nadające się do recyklingu</td></tr> <tr><td>Portuguese</td><td>Embalagem reciclável</td></tr> <tr><td>Romanian</td><td>Ambalaj reciclabil</td></tr> <tr><td>Serbian</td><td>Pakovanje se može reciklirati</td></tr> <tr><td>Slovak</td><td>Recyklovateľné obaly</td></tr> <tr><td>Slovenian</td><td>Embalaža, ki jo je mogoče reciklirati</td></tr> <tr><td>Spanish</td><td>Envase reciclable</td></tr> <tr><td>Swedish</td><td>Återvinningsbar förpackning</td></tr> <tr><td>Turkish</td><td>Geri dönüştürülebilir ambalaj</td></tr> </table>	American English	Recyclable packaging	Bosnian	Pakovanje se može reciklirati	British English	Recyclable packaging	Bulgarian	Рециклируема опаковка	Croatian	Pakiranje koje se može reciklirati	Czech	Recyklovatelný obal	Danish	Genanvendelig emballage	Dutch	Recyclebare verpakking	Estonian	Taaskasutatav pakend	Finnish	Kierrätettävä pakkausmateriaali	French	Emballage recyclable	German	Recyclingfähige Verpackung	Greek	Ανακυκλώσιμη συσκευασία	Hungarian	Újrahasznosítható csomagolás	Icelandic	Endurvinnanlegur umbúðir	Italian	Confezione riciclabile	Latvian	Pārstrādājams iepakojums	Lithuanian	Perdirbama pakuotė	Macedonian	Амбалажата може да се рециклира	Norwegian	Resirkulerbar pakning	Polish	Opakowanie nadające się do recyklingu	Portuguese	Embalagem reciclável	Romanian	Ambalaj reciclabil	Serbian	Pakovanje se može reciklirati	Slovak	Recyklovateľné obaly	Slovenian	Embalaža, ki jo je mogoče reciklirati	Spanish	Envase reciclable	Swedish	Återvinningsbar förpackning	Turkish	Geri dönüştürülebilir ambalaj	<p>To indicate that the marked item or its material is part of a recovery or recycling process</p>	<p>ISO 7000 Graphical symbols for use on equipment — Registered symbols</p>	<p>1135</p>
American English	Recyclable packaging																																																													
Bosnian	Pakovanje se može reciklirati																																																													
British English	Recyclable packaging																																																													
Bulgarian	Рециклируема опаковка																																																													
Croatian	Pakiranje koje se može reciklirati																																																													
Czech	Recyklovatelný obal																																																													
Danish	Genanvendelig emballage																																																													
Dutch	Recyclebare verpakking																																																													
Estonian	Taaskasutatav pakend																																																													
Finnish	Kierrätettävä pakkausmateriaali																																																													
French	Emballage recyclable																																																													
German	Recyclingfähige Verpackung																																																													
Greek	Ανακυκλώσιμη συσκευασία																																																													
Hungarian	Újrahasznosítható csomagolás																																																													
Icelandic	Endurvinnanlegur umbúðir																																																													
Italian	Confezione riciclabile																																																													
Latvian	Pārstrādājams iepakojums																																																													
Lithuanian	Perdirbama pakuotė																																																													
Macedonian	Амбалажата може да се рециклира																																																													
Norwegian	Resirkulerbar pakning																																																													
Polish	Opakowanie nadające się do recyklingu																																																													
Portuguese	Embalagem reciclável																																																													
Romanian	Ambalaj reciclabil																																																													
Serbian	Pakovanje se može reciklirati																																																													
Slovak	Recyklovateľné obaly																																																													
Slovenian	Embalaža, ki jo je mogoče reciklirati																																																													
Spanish	Envase reciclable																																																													
Swedish	Återvinningsbar förpackning																																																													
Turkish	Geri dönüştürülebilir ambalaj																																																													